




TECHNICAL DATA SHEET

DEEP WELL PLATE







MID-KF05US-9-N

Sterile available upon request

GENERAL SPECIFICATIONS					
NUMBER OF WELLS:	96	96			
WELL TOP SHAPE:	Square	Square			
WELL BOTTOM SHAPE:	Round - V	Round - V			
NOTCH LOCATION:	A1	A1			
FOOTPRINT:	ANSI/SLAS 1-2004	ANSI/SLAS 1-2004			
MATERIAL:	Polypropylene	Polypropylene			
COLOR:	Natural	Natural			
TEMPERATURE RANGE:	-80°C to 120°C	-80°C to 120°C			
CENTIFUGATION (MAX.):	4,000 x g	4,000 x g			
AUTOCLAVABLE:	121C at 15 psi, 15 min.	Pre-sterilized			
AUTOMATION COMPATIBLE:	Yes	Yes			

CRITICAL DIMENSIONS				
PLATE HEIGHT:	15.10mm	15.10mm		
PLATE BOTTOM (L x W):	127.50mmx 85.35mm	127.50mmx 85.35mm		
PLATE TOP (L x W):				
FLANGE HEIGHT:	2.50mm	2.50mm		
WELL A1 (X,Y):	14.38mm, 11.24mm	14.38mm, 11.24mm		
WELL DEPTH:	13.30mm	13.30mm		
WELL DIA. (TOP):	8.20mm x 8.20mm	8.20mm x 8.20mm		
WELL DIA. (BOT):				
WELL BOTTOM TO DECK:	1.80mm	1.80mm		
OVERALL VOLUME:	350ul	350ul		
WORKING VOLUME:	250ul	250ul		
VOLUME WITH MAT:				

PACKAGING CONFIGURATIONS PACKAGING QUANTITY (CS): 50; 5 per Pack 50; 5 per Pack CASE DIMENSIONS (IN.): 6 x 10.75 x 17.75 6 x 10.75 x 17.75 CASE WEIGHT (LBS.): 4.5 4.5

STORAGE

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone or ignition.







COMPATIBLE SEALING MAT OPTIONS

ITEM NUMBER:	DESCRIPTION:
MID-MATS-9-PRT	Silicone 96 well, Square Printed
MID-MATS-9-BULK	Silicone 96 well, Square
MID-SQ-9-PRT	Pre-slit Silicone 96, Square Printed

All plates can be sealed with adhesive or heat seal films.