


LF141WWW/OM CUTSHEET




115VAC/60Hz/1PH
NEMA-5-15P

Interior LF141WWW/OM



Exterior LF141WWW/OM



NORLAKE
SCIENTIFIC®

Cu. Ft	Defrost	Door	Int Door	Shelves	Interior Dim.			Refrigerant	H.P.	Amps	Weight
					W"	D"	H"				
14	Manual	1 Solid	0	3 Fixed	29 7/8	32 3/4	62	R134a	1/3	5	125 lbs.

Note: Interior dimensions are estimates and subject to change without notice