


LF161WWW/0M CUTSHEET




115VAC/60Hz/1PH
NEMA-5-15P

Interior LF161WWW/0M



Exterior LF161WWW/0M



NORLAKE
SCIENTIFIC®

Cu. Ft	Defrost	Door	Int Door	Shelves	Interior Dim.			Refrigerant	H.P.	Amps	Weight
					W"	D"	H"				
17	Manual	1 Solid	0	3 Fixed	34	30 3/4	67 1/2	R134a	1/3	5	167 lbs.

Note: Interior dimensions are estimates and subject to change without notice