

Special Offers for Selected Labnet Products

LabValues 2023

From August 1 to December 31, 2023



LV-1.23

Labnet MultiGene™ OptiMax

The MultiGene OptiMax Thermal Cycler delivers advanced speed and features while providing premium performance at an affordable price.

- Unique flexible programming with FlexTemp™ technology
- Fast ramp rates 5°C heating / 3.5°C cooling
- Simple user interface

Cat. No.	Description	List Price	Special Price
TC9610	MultiGene OptiMax with 96-well block, 120V	\$4,118.00	\$3,088.50

Labnet Enduro™ Power Supply

The Enduro Power Supplies are ideal for DNA, RNA and protein electrophoresis. They provide a high degree of sophistication and versatility yet are extremely easy to use and simple to program.

- Dependable performance, intelligent control
- Compact, lightweight and stackable

Cat. No.	Description	List Price	Special Price
E0303	Enduro 300V power supply	\$613.58	\$460.19
E0203	Enduro 250V power supply	\$1,184.14	\$888.11



LV-2.23



LV-3.23

Labnet Enduro UV Transilluminators

The unique hinge design and frame on the UV shield prevents the lid from cracking and hinge damage predominant on most transilluminators in the market.


- Compact size with 21 x 26 cm viewing area
- Available in single or dual wavelength models

Cat. No.	Description	List Price	Special Price
U1001	Enduro UV Transilluminator with 302nm wavelength, 115V	\$1,880.49	\$1,128.29
U1002	Enduro UV Transilluminator with 302nm and 365 nm wavelengths, 115V	\$1,985.52	\$1,290.59

Buy Labnet AccuBlock™ Digital Dry Bath and Get Block for FREE



LV-4.23



Purchase a Labnet AccuBlock Digital Dry Bath (**single block**) and **get one (1) block at no additional cost.**

Purchase a Labnet AccuBlock Digital Dry Bath (**dual block**) and **get two (2) blocks at no additional cost.**

Cat. No.	Description
D1301	Digital dry bath, single block, 120V
D1302	Digital dry bath, dual block, 120V

Get One or Two Blocks for FREE

Cat. No.	Description
D1101	Solid block for machining (no holes)
D1102	Single block, 48 x 0.2 mL PCR tubes or 6 x 0.2 mL
D1102A	Single block, 24 x 2.0 mL tubes
D1105	Single block, 24 x 0.5 mL tubes
D1105A	Single block, 24 x 1.5 mL tubes
D1106	Single block, 35 x 6 mm tubes
D1110	Single block, 20 x 10 mm tubes
D1112	Single block, 20 x 12 mm tubes

Cat. No.	Description
D1113	Single block, 20 x 13 mm tubes
D1115-TALL	Single block, 12 x 15 mL centrifuge tubes
D1116	Single block, 12 x 15 or 16 mm tubes
D1117	Single block, 12 x 17 mm tubes
D1120	Single block, 6 x 20 mm tubes
D1125	Single block, 6 x 25 mm tubes
D1150-TALL	Single block, 5 x 50 mL conical tubes
D1196-PCR	Single block, 96-well PCR plate, skirted or non-skirted
D1296	Dual block, 96-well microplate or 4 slides
D1296-PCR	Dual block, 96-well PCR plate, skirted or non-skirted
D12384	Dual block, 384-well PCR plate (for dual block unit only)

—\$580

LV-5.23



Labnet Benchtop Shaking Incubator 222DS

The 222DS lab incubator lid has a low pivot point which provides a compact lab friendly shaking incubator sure to fit on the lab bench and under your shelves.

- ◉ Exceptional versatility
- ◉ Precise SmartCheck™ digital microprocessor temperature control

Cat. No.	Description	List Price	Special Price
I-5222-DS	222DS Benchtop Shaking Incubator, 120V	\$3,867.51	\$3,287.38


Labnet Enduro™ GDS II Gel Documentation Systems

Enduro GDS II and GDS Touch II are used for visualization of DNA and protein stained within a gel and provide photo documentation of these gels in a laboratory setting.

The GDS II comes in either a 302 nm or 365 nm version and runs on a Windows-based PC (not included). The GDS Touch II comes in either a 302 nm or 365 nm wavelength version with 470 nm epi-blue lights for use with "safe dyes" that eliminate the need for ethidium bromide.

⦿ Region of interest function allows optimal exposure in less than 30 seconds. No more multiple exposures to get good images.

Cat. No.	Description	List Price	Special Price
GDS2-1302	Enduro GDS II with 302 nm UV	\$6,073.92	\$4,859.14
GDS2-1365	Enduro GDS II with 365 nm UV	\$6,073.92	\$4,859.14
GDST2-1302	Enduro GDS Touch II with epi-blue and 302 nm UV	\$9,874.08	\$7,899.26



Save BIG with Labnet LabValues Special Offers



LV-7.23

Labnet Excel™ Electronic Pipettors

The Excel is an ergonomic, intuitive electronic pipettor for use in most laboratory applications. Its high precision stepping motor controls piston action, resulting in exceptional accuracy and precision.

- ⦿ Factory calibrated to strict ISO8655 standard
- ⦿ Large LCD screen and simple interface
- ⦿ Ergonomically designed handle and finger rest

Cat. No.	Description	List Price	Special Price
P3600L-200	20-200 µL Excel electronic pipettor, 1-channel	\$394.29	\$315.43
P3600L-1200	100-1,200 µL Excel electronic pipettor, 1-channel	\$394.29	\$315.43
P3608L-10	0.5-10 µL Excel electronic pipettor, 8-channels	\$805.12	\$644.10
P3608L-20	2-20 µL Excel electronic pipettor, 8-channels	\$805.12	\$644.10
P3608L-200	20-200 µL Excel electronic pipettor, 8-channels	\$805.12	\$644.10
P3608L-1200	100-1,200 µL Excel electronic pipettor, 8-channels	\$912.70	\$730.16
P3612L-10	0.5-10 µL Excel electronic pipettor, 12-channels	\$939.79	\$751.83
P3612L-20	2-20 µL Excel electronic pipettor, 12-channels	\$939.79	\$751.83
P3612L-200	20-200 µL Excel electronic pipettor, 12-channels	\$939.79	\$751.83
P3612L-1200	100-1,200 µL Excel electronic pipettor, 12-channel	\$1073.88	\$859.10

Labnet BioPette™ Plus Pipettors

BioPette Plus Pipettors feature smooth plunger movement with low pipetting forces. Volume can be set by turning either the thumb wheel or the plunger button and is protected by a locking mechanism.

All BioPette Plus Pipettors feature an accurate and precise pipetting mechanism, enclosed in an ergonomic, well-balanced, and robust handle.

- Factory calibrated to strict ISO 8655 standards
- Calibration key is included with each pipettor for easy in-house calibration
- Fully autoclavable and UV-resistant
- Designed for universal fit with all common brands of pipet tips
- 3-year warranty

Cat. No.	Description	List Price	Special Price
P3942-10	0.5-10 µL Biopette Plus pipettor, single-channel	\$273.87	\$219.10
P3942-20	2-20 µL Biopette Plus pipettor, single-channel	\$257.90	\$206.32
P3942-100	10-100 µL Biopette Plus pipettor, single-channel	\$257.90	\$206.32
P3942-200	20-200 µL Biopette Plus pipettor, single-channel	\$257.90	\$206.32
P3942-1000	100-1,000 µL Biopette Plus pipettor, single-channel	\$257.90	\$206.32
P3942-SK4	Labnet BioPette Plus starter kit	\$933.91	\$747.13



Labnet Labpette™ FX Fixed Volume Pipettors

Labnet Labpette™ FX Fixed Volume Pipettors are ideal for applications where specific volumes are used repeatedly and sterilization by autoclaving is required.

- Factory calibrated to strict ISO 8655 standards
- Seventeen fixed volumes from 2 to 1,000 µL
- Ergonomically designed handle and finger rest

Cat. No.	Description	List Price	Special Price
P3900-2	2 µL Labpette FX fixed volume pipettor	\$114.10	\$91.28
P3900-5	5 µL Labpette FX fixed volume pipettor	\$112.97	\$90.38
P3900-10	10 µL Labpette FX fixed volume pipettor	\$114.10	\$91.28
P3900-50	50 µL Labpette FX fixed volume pipettor	\$115.23	\$92.18
P3900-200	200 µL Labpette FX fixed volume pipettor	\$116.38	\$93.10
P3900-250	250 µL Labpette FX fixed volume pipettor	\$115.23	\$92.18
P3900-300	300 µL Labpette FX fixed volume pipettor	\$115.23	\$92.18
P3900-400	400 µL Labpette FX fixed volume pipettor	\$112.97	\$90.38
P3900-500	500 µL Labpette FX fixed volume pipettor	\$116.38	\$93.10



-33%

Laboratory Plastics and Consumables

Cat. No.	Description	List Price	Special Price
K1941	Semi-micro UV range cuvette	\$118.02	\$79.07
K1938	Semi-micro VIS range cuvette	\$130.25	\$87.27
K1962	Cuvette caps, 1000/pk	\$58.20	\$38.99
K1961	UV range cuvette with 4 clear sides, PMMA, 4.5 mL	\$147.11	\$98.56
K1960	VIS range cuvette with 4 clear sides, PS, 4.5 mL	\$106.14	\$71.11
K1720	Graduated beaker PMP 25 mL, 40/pk	\$207.30	\$138.89
K1721	Graduated beaker PMP 50 mL, 40/pk	\$270.46	\$181.21
K1723	Graduated beaker PMP 250 mL, 16/pk	\$149.74	\$100.33
K1160	Graduated beaker with handle short, PP, 5000 mL, 6/pk	\$311.02	\$208.38
K1415	Untapered 4-way connector, 13.5 mm	\$353.56	\$236.89
K1253	Untapered L connector, 10 mm	\$130.25	\$87.27
P9615-50	Microplate, 96 x 1.0 mL/round well, polypropylene, nonsterile	\$328.34	\$219.99
P9616-50	Microplate, 96 x 1.0 mL/round well, polypropylene, sterile	\$363.64	\$243.64
P9635	Microplate, 96 x 2.0 mL/square well, polypropylene, nonsterile	\$354.96	\$237.82
P9639	Sealing mat for 96 square well microplate, sterile	\$212.27	\$142.22
P9651	384-well microplate, 384 x 240 µL, sterile	\$425.96	\$285.39
P9702	High volume deep well microplates, 48 x 5 mL well, SBS format, 44 mm high	\$371.56	\$248.95
R8000	Styrofoam storage box, 100 x 1.5/2.0 mL	\$8.62	\$5.78
R8315	Vertical freezer rack for 9 boxes, fits chest type freezers (150 x 150 x 510 mm)	\$197.68	\$132.45
R8300	100-place microtube box with hinged lid, natural	\$10.50	\$7.04
R8300-Y	100-place microtube box with hinged lid, yellow	\$11.33	\$7.59
R8300-B	100-place microtube box with hinged lid, blue	\$11.33	\$7.59
R8300-O	100-place microtube box with hinged lid, orange	\$11.33	\$7.59
R8300-G	100-place microtube box with hinged lid, green	\$11.33	\$7.59
R8300-BLK	100-place microtube box with hinged lid, opaque black for light sensitive samples	\$12.15	\$8.14

**SAVE
BIG**

**Special offers valid
August 1 – December 31, 2023**