

# Leave the Ice Buckets Behind!

## IsoFreeze® Racks are BETTER than ice.



Maintain sample temp less than 4°C for up to 4 hour when used at room temperature.

- Contrasting colors indicate safe sample temperatures at a glance.
- Patented color-changing material indicates when samples are approaching 7°C.

TEMPERATURE BELOW 4°C

TEMPERATURE CHANGING

TEMPERATURE AT 7°C



### MCT Color-Change Racks

- Holds 24 x 0.5, 1.5 or 2.0 mL tubes.

5630-T4	- Cold Temperature, Purple - Warm Temperature, Pink	146 x 105 x 70 mm 5.75 x 4.13 x 2.76 in.	1/pk	\$79.78
5630-T6	- Cold Temperature, Green - Warm Temperature, Yellow	146 x 105 x 70 mm 5.75 x 4.13 x 2.76 in.	1/pk	\$79.78

### PCR Color-Change Racks

- Holds 96 x 0.2 or 0.1 mL tubes, strip tubes, or 96 wells plates.

5640-T4	- Cold Temperature, Purple - Warm Temperature, Pink	141 x 99 x 44 mm 5.55 x 3.90 x 1.73 in.	2/pk	\$93.43
5640-T6	- Cold Temperature, Green - Warm Temperature, Yellow	141 x 99 x 44 mm 5.55 x 3.90 x 1.73 in.	2/pk	\$93.43

## IsoFreeze® Flipper® Racks

Maintain cold temperature for hours on the bench!

Side A: holds 20 x 2.0 mL tubes

Side B: holds 20 x 0.5 mL tubes



### Blue -20°C IsoFreeze Flipper Racks

- Maintain temp between -10°C and -20°C for up to 3 hours.

TFL-20	Blue, -20°C IsoFreeze Flipper Rack	1/pk	\$104.22
--------	------------------------------------	------	----------

### White 0°C IsoFreeze Flipper Racks

- Maintain temp of 0°C for up to 5 hours.

TFL-0	White, 0°C IsoFreeze Flipper Rack	1/pk	\$104.22
-------	-----------------------------------	------	----------

## New! Microcentrifuge Storage Racks

MIDSCI is pleased to introduce a new 5 mL tube rack, expanding our line of handy storage solutions for most tubes in the lab.

- Three models accommodate common tube sizes (0.5, 1.5/2.0, and 5.0 mL).
- Ideal for screw cap microtubes and microcentrifuge tubes.
- Small footprint saves bench and freezer space.
- Side notches facilitate easy lid removal.
- Suitable for temperatures as low as -90°C as well as autoclaving.
- Racks come with lid.



5527-00	5 mL 21-Place Freezer Box, Natural 143.8 x 75.2 x 53.0 mm   5.66 x 2.96 x 2.09 in.	20/cs	\$105.78
5520-00	1.5 - 2.0 mL 50-Place Freezer Box, Natural 143.8 x 75.2 x 53.0 mm   5.66 x 2.96 x 2.09 in.	20/cs	\$118.66
5530-00	0.5 mL 50-Place Freezer Box, Natural 139.7 x 71.1 x 61.8 mm   5.5 x 2.80 x 2.43 in.	20/cs	\$118.66

# Check Out Our Best Selling Racks!

## Freezer Racks

Fits in standard freezer racks!\*

### 81-Place Freezer Racks

- Suitable for temperatures as low as -90°C
- Sturdy polypropylene
- Stores tubes up to 2.0 mL
- Bright colors facilitate organization



### 100-Place Freezer Racks

- Hinged lid
- Holds 1.5/2.0 mL tubes & 1.2 mL cryo vials
- Rated to -90°C
- 5.55 x 6.05 x 2.15 in.



FB-81*	81-Place Freezer Rack w/ Lift-off Lid, Natural	5/pk	\$36.91
FB-81C*	81-Place Freezer Rack w/ Lift-off Lid, Assorted Fluorescent Colors	5/pk	\$79.78

FRB-100	100-Place Freezer Rack, Natural	5/pk	\$41.14
FRB-100C	100-Place Freezer Rack, Assorted Fluorescent	5/pk	\$41.14

Contact us for individual colors!

Contact us for individual colors!

### CryoFreeze® Cryo Storage Racks

- Durable polycarbonate
- Drains & vents prevent pressure build-up
- Crystal clear lid
- Easy-to-read alphanumeric grid



### 15 mL & 50 mL Freezer Racks

- Polypropylene
- Removable base for use in ice bath
- 5.38 x 5.38 x 5 in. with lid




CF-25*	25-Place Rack, Holds 1.0/1.8 mL, Assorted Colors	5/pk	\$39.16
CF-81*	81-Place Rack, Holds 1.0/1.8 mL, Assorted Colors	5/pk	\$59.97
CF-81-LG	81-Place Rack, Holds 3.6/4.5 mL, Assorted Colors	5/pk	\$80.69
CF-100*	100-Place Rack, Holds 1.0/1.8 mL, Assorted Colors	5/pk	\$65.03

5570-02	36-Place Freezer Rack, Holds 36 x 15 mL, Green	2/pk	\$37.44
5580-01	16-Place Freezer Rack, Holds 16 x 50 mL, Blue	2/pk	\$37.44

## Utility Racks

### 96-Place PCR Work Up Racks


- Working storage rack for PCR processing
- Standard spacing for strip tubes or plates



WRI-4A



WR-96A




WUR-96A

WRI-4A	PCR WorkUp Rack & Workstation w/ Lid, Assorted Fluorescent	5/pk	\$126.81
WR-96A	96-Place PCR WorkUp Rack, Assorted Fluorescent	5/pk	\$39.27
WUR-96A	96-Place Stacking PCR WorkUp Rack, Assorted Standard	5/pk	\$42.36

### 80-Place Storage Racks

- Versatile & efficient storage
- Fraction collector plate style
- Suitable for tubes up to 2.0 mL



HS29025A	80-Place Storage Rack, Natural	5/pk	\$17.68
HS29025G	80-Place Storage Rack, Assorted Fluorescent	5/pk	\$19.80

### 96-Place Flipper® Racks for PCR


- One side holds 0.5 mL and the other 1.5/2.0 mL tubes
- Bright colors enable organization
- 8 x 12 grid
- Also available in 8 x 8 and 4 x 8 formats



FL-96F-5R	96-Place Flipper Rack, Clear Lid, Assorted Fluorescent	5/pk	\$51.09
FL-96CL-B-5R	96-Place Flipper Rack, Clear Lid, Blue Standard	5/pk	\$60.22
FL-96C-5R	96-Place Flipper Rack w/ Clear Lid, Assorted Standard	5/pk	\$51.09
FL-96RO-A-5R	96-Place Flipper Rack only, Assorted Standard	5/pk	\$44.22

### 4-Way Flipper® Racks

- Rugged polypropylene
- Bright colors
- Autoclavable
- Suitable for ULT freezers (to -90°C)



5400-29	4-Way Flipper, Interlocking, Assorted Fluorescent. Holds: 32 x 0.5 mL, 32 x 1.5/2.0 mL, 12 x 15, 4 x 50 mL tubes	4 rks/tray	\$59.90
5400-39	4-Way Flipper, Assorted Standard Colors. Holds: 40 x 1.5/2.0 mL, 24 x 15 mL, 20 x 20 mm, 12 x 25 mm Tubes	4 rks/tray	\$56.84