

**SPEC SHEET ABT-30R**

<b>Product Name:</b>	ABT-30R
<b>Features:</b>	<ul style="list-style-type: none"> <li>• General Purpose Refrigerator</li> <li>• 30 Cu. Ft. Capacity</li> <li>• One solid swing door, right hinged, not reversible</li> <li>• Mechanical thermostat</li> <li>• Temperature range: 0°C to 10°C</li> <li>• Five adjustable shelves</li> <li>• Keyed door lock</li> <li>• Magnetic door gasket for positive seal</li> <li>• Stainless steel door handle</li> <li>• Leveling legs</li> <li>• Probe access port (3/8")</li> <li>• Powder coat paint exterior</li> <li>• High density urethane foam cabinet and door insulation</li> <li>• CFC-free refrigerant (R134a), foam insulation and packaging</li> <li>• Two year parts and labor warranty, plus an additional three year compressor parts warranty</li> <li>• Exterior dimensions: 35 ¾"W x 35" D x 72 ½"H</li> <li>• Auto defrost</li> <li>• 115V, 60Hz, 6.8 amps, 1/3 HP</li> <li>• UL/C-UL listed</li> <li>• Shipping weight: 396 lbs.</li> </ul>