

SPEC SHEET ABT-HC-10PG

Product Name:	ABT-HC-10PG
Features:	<ul style="list-style-type: none"> • Premier Glass Door Compact Laboratory Refrigerator • 10.5 Cu. Ft. Capacity • One glass swing door, right hinged, not reversible • Microprocessor temperature controller • Digital temperature display with sample simulated probe • Audible and visual high and low temperature alarms • Remote alarm contacts • Adjustable operating temperature range: 1°C to 10°C • Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall • Bottom footed shelf to allow air circulation • Forced draft air circulation • Back wall evaporator cover • Keyed door lock • Magnetic door gasket for positive seal • Leveling legs • LED interior light • Probe access port (3/8") • Powder coat exterior • High density urethane foam cabinet insulation • HFC-free refrigerant (R600a), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts warranty • Exterior dimensions: 23 ¾"W x 26 ½"D x 59"H • Note: This unit must have 4" clearance on sides and back for adequate ventilation • UL/C-UL listed • ENERGY STAR CERTIFIED • Cycle defrost • 115V, 60 Hz, 1/3 HP, 2.3 Amps • Energy Consumption 2.42 kWh/day • Shipping weight: 209 lbs.