

SPEC SHEET ABT-HC-19-LH

Product Name:	ABT-HC-19-LH
Features:	<ul style="list-style-type: none"> • Premier Glass Door Laboratory Refrigerator • 19 Cu. Ft. capacity • One swing glass door, self-closing (energy efficient, double pane vacuum insulated), LEFT hinged, not reversible • Microprocessor Temperature Controller • Adjustable operating temperature range: 1°-10°C • Standard Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms • Sample simulated digital temperature display; .1° Resolution • Adjustable temperature display in Celsius or Fahrenheit • Audible and visual high/low temperature alarms • Sensor error alarm, Power failure alarm • Resettable 1-point Min/Max History • Remote alarm contacts • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • Four adjustable shelves (adjustable in 1/2" increments) • Forced draft air circulation • Keyed door lock, Magnetic door gasket for positive seal, Recessed door handle • LED interior lights are safety shielded and switch controlled • Swivel casters (adds 4" to height), Probe access port (3/8") • White laminated vinyl exterior, White aluminum liner and stainless steel floor interior • HFC-free refrigerant (R290), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts warranty • Exterior dimensions: 27 1/8" W x 25 1/2" D x 82 5/8" H with casters; 78 5/8" H without casters • Cycle defrost • 115V, 60 Hz, 5.4 Amps, 1/3 HP, BTU/HR 4932 • UL/C-UL listed • ENERGY STAR CERTIFIED • Shipping weight: 295 lbs.