

SPEC SHEET **ABT-HC-LP-33-TS**

Product Name:	ABT-HC-LP-33-TS
Features:	<ul style="list-style-type: none"> • Templog Premier Glass Door Laboratory Refrigerator • 33 Cu. Ft. capacity • Two sliding glass doors, self-closing (energy efficient, double pane vacuum insulated) • Microprocessor temperature controller • Adjustable operating temperature range: 1°C - 10°C • Touch Screen Digital Temperature Display and Alarm Module • Password-protected control panel • Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit • Active temperature scale bar display • Resettable 1-point min/max history • Digital 24-hour temperature chart and Temperature data logging • USB port for data transfer and software updates • Real-time operational display clock • 4-20mA output • Electronic troubleshooting guide and Electronic QR Code to access Owner's Manual & technical support contact information • Adjustable alarm volume • Audible and visual high/low temperature alarms 2 probes; 1 in air and 1 in sample bottle • Sensor error alarm; Power failure alarm; Door ajar alarm • Low battery alarm • Operation indicator • Alarm validation (audible, visual, and remote alarm) and alarm mute with ring back

SPEC SHEET ABT-HC-LP-33-TS

Product Name:	ABT-HC-LP-33-TS
Features:	<ul style="list-style-type: none"> • Remote alarm contacts • Eight adjustable shelves (adjustable in ½" increments) • Forced draft air circulation • Keyed door locks • Magnetic door gasket for positive seal • Integrated door handles • LED interior lights are safety shielded and switch controlled • Swivel Casters (3 ½"), Probe access port (¾") • White powder-coated steel exterior and interior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) • Three year parts and labor warranty, plus an additional two year compressor parts warranty • Exterior dimensions: 39 ⅝" W x 33 ¼" D x 81 ¾" H with casters; 78 ¼" H without casters • Cycle defrost • 115V, 60 Hz, 4.5 Amps, 1/4 HP • ETL listed • Shipping weight: 419 lbs.