

**SPEC SHEET**    **ABT-HC-MFP-17**

<b>Product Name:</b>	ABT-HC-MFP-17
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Standard Natural Refrigerant Manual Defrost Laboratory Freezer</li> <li>• 17 Cu. Ft. Capacity</li> <li>• One solid swing door, right hinged</li> <li>• Microprocessor temperature controller</li> <li>• Digital temperature display</li> <li>• Audible and visual high and low temperature alarms</li> <li>• Remote alarm contacts</li> <li>• 2 probes; 1 in air and 1 in glycol bottle</li> <li>• Adjustable operating temperature range: -15°C – -25°C</li> <li>• Six fixed shelves</li> <li>• Keyed door lock</li> <li>• Magnetic door gasket for positive seal</li> <li>• Leveling legs</li> <li>• Probe access port (3/8")</li> <li>• Powder coat exterior</li> <li>• High density urethane foam cabinet and door insulation</li> <li>• Natural, Hydrocarbon Refrigerant (R290)</li> <li>• Two year parts and labor warranty, plus additional three year compressor parts warranty</li> <li>• Exterior dimensions: 30 3/4" W x 28 5/8" D x 67 3/4" H</li> <li>• Note: This unit must have 4" clearance on sides and back for adequate ventilation</li> <li>• Manual defrost</li> <li>• 115V, 60 Hz, 6 Amps, 1/2 HP</li> <li>• ETL listed</li> <li>• Shipping weight: 220 lbs.</li> </ul>