

**SPEC SHEET**    **ABT-HC-RFC7S**

<b>Product Name:</b>	ABT-HC-RFC7S
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Refrigerator &amp; Freezer Combination</li> <li>• 5.2 cu. ft. capacity (Refrigerator) / 1.7 cu. ft. capacity (Freezer)</li> <li>• Two exterior swing doors (solid door refrigerator/solid door freezer), right hinged, not reversible</li> <li>• Microprocessor temperature controller, refrigerator and freezer</li> <li>• Digital temperature display, refrigerator and freezer</li> <li>• Audible and visual high and low temperature alarms, refrigerator and freezer</li> <li>• 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer</li> <li>• Remote alarm contacts, refrigerator and freezer</li> <li>• Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer:-15°C – -25°C</li> <li>• 3 adjustable shelves (refrigerator) / 2 adjustable shelves (freezer), shelves have back rails to prevent product from being pushed against back wall</li> <li>• Footed bottom shelf refrigerator, and freezer</li> <li>• Keyed door locks</li> <li>• Magnetic door gasket for positive seal</li> <li>• Integrated door handles</li> <li>• Leveling legs Stacking kit</li> <li>• Two probe access ports (3/8"); one, refrigerator/one, freezer</li> <li>• White powder coated exterior</li> <li>• HFC-free refrigerant (R600a), foam insulation and packaging</li> <li>• Two year parts and labor warranty, plus an additional three year compressor parts warranty</li> </ul>

**SPEC SHEET**    **ABT-HC-RFC7S**

<b>Product Name:</b>	ABT-HC-RFC7S
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Exterior dimensions: 23 3/4" W x 24" D x 53 1/8" H</li> <li>• Manual defrost freezer / cycle defrost refrigerator</li> <li>• 115V, 60 Hz, 1.3 Amps refrigerator, 1 Amp freezer, 1/6 HP refrigerator, 1/8 HP freezer</li> <li>• UL/C-UL listed</li> <li>• Shipping weight: 215 lbs.</li> <li>• *Please note unit requires 2 outlets.</li> </ul>