

**SPEC SHEET**    **ABT-HC-RFC9G**

<b>Product Name:</b>	ABT-HC-RFC9G
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Refrigerator &amp; Freezer Combination</li> <li>• 9 Cu. Ft. Capacity (4.6 Ref./4.2 Freezer)</li> <li>• Two exterior swing doors (glass door refrigerator/solid door freezer), right hinged*</li> <li>• Microprocessor temperature controller, refrigerator and freezer</li> <li>• Digital temperature display, refrigerator and freezer</li> <li>• Audible and visual high and low temperature alarms, refrigerator and freezer</li> <li>• 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer</li> <li>• Remote alarm contacts, refrigerator and freezer</li> <li>• Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer: -15°C - -25°C</li> <li>• 2 adjustable/1 fixed shelves refrigerator with guard rail on back to prevent product from being pushed against back wall/2 fixed shelves freezer</li> <li>• Keyed door locks</li> <li>• Magnetic door gasket for positive seal</li> <li>• Recessed door handle, freezer; heavy duty door handle, refrigerator</li> <li>• Leveling legs</li> <li>• Stacking kit</li> <li>• LED interior light, Refrigerator</li> <li>• Two probe access ports (3/8"); one, refrigerator/one, freezer</li> <li>• White powder coated exterior</li> <li>• High density urethane foam cabinet and door insulation</li> <li>• HFC-free refrigerant (R600a), foam insulation and packaging</li> </ul>

**SPEC SHEET**    **ABT-HC-RFC9G**

<p><b>Product Name:</b></p>	<p>ABT-HC-RFC9G</p>
<p><b>Features:</b></p>	<ul style="list-style-type: none"> <li>• Exterior dimensions: 23 ¾" W x 25 ½" D x 66 ¾" H</li> <li>• Two year parts and labor warranty, plus an additional three year compressor parts warranty</li> <li>• Manual defrost freezer / cycle defrost refrigerator</li> <li>• 115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer</li> <li>• UL/C-UL listed</li> <li>• Shipping weight: 282 lbs.</li> <li>• *Please note this unit requires 2 outlets.</li> <li>• *This product is available left hinged, product number: ABT-HC-RFC9G-LH</li> </ul>