

SPEC SHEET ABT-HC-SLS-16

| | |
|----------------------|---|
| Product Name: | ABT-HC-SLS-16 |
| Features: | <ul style="list-style-type: none"> • Standard Solid Door Laboratory Refrigerator • 16 Cu. Ft. capacity • One swing solid door, self-closing, right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Remote alarm contacts • Adjustable operating temperature range: 1°C – 10°C • Three adjustable shelves (adjustable in ½" increments) • Forced draft air circulation • Heavy duty door handle • Keyed door lock • Magnetic door gasket for positive seal • LED interior lights are safety shielded and switch controlled • Swivel Casters (3 ½"); Probe access port (3/8") • White powder-coated steel exterior and interior • Hydrocarbon, natural refrigerant (R290) • One year parts and labor warranty, plus an additional four year compressor parts warranty • Exterior dimensions: 25" W x 29 ¾" D x 79" H with casters; 75 ½" without casters • Cycle defrost • 115V, 60 Hz, 3 Amps, 1/5 HP • ETL listed • Shipping weight: 277 lbs. |