

**SPEC SHEET ABT-HC-SLS-26**

<b>Product Name:</b>	ABT-HC-SLS-26
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Standard Solid Door Laboratory Refrigerator</li> <li>• 26 Cu. Ft. capacity</li> <li>• One swing solid door, self-closing, right hinged, not reversible</li> <li>• Microprocessor temperature controller</li> <li>• Digital temperature display</li> <li>• Audible and visual high and low temperature alarms</li> <li>• Remote alarm contacts</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Adjustable operating temperature range: 1°C – 10°C</li> <li>• Four adjustable shelves (adjustable in ½" increments)</li> <li>• Forced draft air circulation</li> <li>• Keyed door lock</li> <li>• Magnetic door gasket for positive seal</li> <li>• Heavy duty door handle</li> <li>• LED interior lights are safety shielded and switch controlled</li> <li>• Swivel Casters (3 ½"); Probe access port (¾")</li> <li>• White powder-coated steel exterior and interior</li> <li>• Hydrocarbon, natural refrigerant (R290)</li> <li>• One year parts and labor warranty, plus an additional four year compressor parts warranty</li> <li>• Exterior dimensions: 28 ⅜" W x 36 ½" D x 81 ¾" H with casters; 78 ¼" H without casters</li> <li>• Cycle defrost</li> <li>• 115V, 60 Hz, 3 Amps, 1/5 HP</li> <li>• ETL listed</li> <li>• Shipping weight: 275 lbs.</li> </ul>