

SPEC SHEET ABT-HC-SLS-36

Product Name:	ABT-HC-SLS-36
Features:	<ul style="list-style-type: none"> • Standard Solid Door Laboratory Refrigerator • 36 Cu. Ft. capacity • Two swing solid doors, self-closing, right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • 2 probes; 1 in air and 1 in sample bottle • Adjustable operating temperature range: 1°C – 10°C • Keyed door locks • Magnetic door gasket for positive seal • Eight adjustable shelves (adjustable in ½" increments) • Forced draft air circulation • Heavy duty door handles • LED interior lights are safety shielded and switch controlled • Swivel Casters (3 ½"); Probe access port (⅜") • White powder-coated steel exterior and interior • Hydrocarbon, natural refrigerant (R290) • One year parts and labor warranty, plus an additional four year compressor parts warranty • Exterior dimensions: 39 ⅝" W x 34 ¾" D x 81 ¾" H with casters; 78 ¼" H without casters • Cycle defrost • 115V, 60 Hz, 4.5 Amps, 1/4 HP • ETL listed • Shipping weight: 360 lbs.