

SPEC SHEET ABT-HC-SSP-23G

Product Name:	ABT-HC-SSP-23G
Features:	<ul style="list-style-type: none"> • Premier Series Stainless Steel Glass Door Refrigerator • 23 Cu. Ft. Capacity • One swing glass door, self-closing, right hinged, not reversible • 4 ea. Adjustable Shelves • Digital display microprocessor temperature controller • Adjustable operating temperature: 1°C to 10°C (34°F to 50°F) • Audible and visual Hi/Lo temperature alarm • 2 probes; 1 in air and 1 in sample bottle • Sensor error alarm • Continuous product (simulated) temperature display • Remote alarm contacts • Forced draft air circulation • Keyed door lock • Magnetic door gaskets for positive seal • Heavy duty door handle • LED interior lights are safety shielded and switch controlled • Probe access port (3/4") • Stainless steel exterior • White-powder coated aluminum interior with 304 stainless steel floor • Easily cleaned smooth interior • Hydrocarbon natural refrigerant (R290) • 2 year parts and labor warranty, plus an additional 3 year compressor parts warranty • 4 swivel casters (5 3/4")

SPEC SHEET **ABT-HC-SSP-23G**

Product Name:	ABT-HC-SSP-23G
Features:	<ul style="list-style-type: none">• Exterior dimensions: 26 3/4" W x 35 3/4" D x 83 1/4" H with casters; 77 1/2" H without casters• Cycle defrost• 115V, 60 Hz, 3 Amps, 1/5 HP• ETL listed• Shipping weight: 356 lbs.