

SPEC SHEET **ABT-HC-UCBI-0204-LH**

Product Name:	ABT-HC-UCBI-0204-LH
Features:	<ul style="list-style-type: none"> • Premier Undercounter Refrigerator Built-In • 2.5 Cu. Ft. Capacity • One solid swing door, left hinged • Microprocessor temperature controller • Adjustable operating temperature range: 1°C – 10°C • Digital temperature display • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Remote alarm contacts • Four adjustable shelves (two large, two small) with guard rail on back to prevent product from being pushed against the back wall • Back wall evaporator cover • Forced draft circulation • Self-closing door • Keyed door lock • Magnetic door gasket for positive seal • Leveling legs • Probe access port (3/8") • Powder coat exterior • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R600a), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts only warranty • Pyxis®, Omnicell® and AcuDose RX® compatible • Exterior dimensions: 17 ¾ "W x 22 ¾ "D x 30 ¾"H • Energy Star Certified for Laboratory Refrigerators and Freezers

SPEC SHEET **ABT-HC-UCBI-0204-LH**

Product Name:	ABT-HC-UCBI-0204-LH
Features:	<ul style="list-style-type: none">• Cycle defrost• 115V, 60 Hz, .90 Amps, 1/8 HP• Energy Consumption 0.58 kWh/day• ADA compliant• UL/C-UL listed• Shipping weight: 85 lbs.