

SPEC SHEET ABT-HC-UCFS-0220M

Product Name:	ABT-HC-UCFS-0220M
Features:	<ul style="list-style-type: none"> • Standard Under Counter Freezer Freestanding • 1.5 Cu. Ft. Capacity • One solid swing door, right hinged • Mechanical thermostat • Average cycle temperature: -20°C • One fixed shelf • Keyed door lock • Magnetic door gasket for positive seal • Heavy duty door handle • Leveling legs • Probe access port (3/8") • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R600a), foam insulation and packaging • One year parts and labor warranty • Pyxis®, Omnicell® and AcuDose RX® compatible • Exterior dimensions: 18 ½" W x 19 ½" D x 19 ½" H • Manual defrost • 115V, 60 Hz, 1/10 HP, 1.1 amps • Shipping weight: 51 lbs.