

SPEC SHEET ABT-HC-UCFS-0440

<p>Product Name:</p>	<p>ABT-HC-UCFS-0440</p>
<p>Features:</p>	<ul style="list-style-type: none"> • Premier Undercounter Freezer Freestanding (-40°C) • 4 Cu. Ft. Capacity • One solid swing door, right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Remote alarm contacts • Adjustable operating temperature range: -37°C – -43°C • Three adjustable shelves with guard rail on back to prevent product from being pushed against the back wall • Bottom footed shelf to allow air circulation • Back wall evaporator cover • Keyed door lock • Magnetic door gasket for positive seal • Recessed door handle • Leveling legs • LED interior light • Probe access port (3/8") • White cabinet and door • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R290), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts only warranty • Pyxis®, Omnicell® and AcuDose RX® compatible • Exterior dimensions: 23 ¾" W x 24" D x 32 1/8" H

SPEC SHEET **ABT-HC-UCFS-0440**

Product Name:	ABT-HC-UCFS-0440
Features:	<ul style="list-style-type: none">• Note: This unit must have 4" clearance on sides and back for adequate ventilation• Manual defrost• 115V, 60 Hz, 2.2 Amps, 1/4 HP• Energy Consumption 1.34 kWh/day• UL/C-UL listed• ENERGY STAR CERTIFIED• Shipping weight: 128 lbs.