

SPEC SHEET ABT-HCPTP-49

<p>Product Name:</p>	<p>ABT-HCPTP-49</p>
<p>Features:</p>	<ul style="list-style-type: none"> • Premier Pass Thru Laboratory Refrigerator • 49 Cu. Ft. capacity • Two glass swing doors on front, two solid swing doors on back • Microprocessor Temperature Controller • Adjustable operating temperature range: 1°-10°C • Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms • Sample simulated digital temperature display; .1° Resolution • Adjustable temperature display in Celsius or Fahrenheit • Audible and visual high/low temperature alarms • Sensor error alarm • Power failure alarm • Resettable 1-point Min/Max History • Remote alarm contacts • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • Eight adjustable shelves • Forced draft air circulation • Magnetic door gasket for positive seal • Casters (adds 4" to height) • Heavy duty door handle • One probe access port (3/8") • White textured cabinet and door • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R290), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts warranty

SPEC SHEET **ABT-HCPTP-49**

Product Name:	ABT-HCPTP-49
Features:	<ul style="list-style-type: none">• Exterior dimensions: 54 1/8" W x 33 7/8" D x 82 5/8" H with casters; 78 5/8" H without casters• Cycle defrost• 115V, 60 Hz, 8.5 Amps, 1/2 HP, BTU/HR 7315• UL/C-UL listed• Shipping weight: 575 lbs.