

SPEC SHEET **ABT-MFP-3030**

Product Name:	ABT-MFP-3030
Features:	<ul style="list-style-type: none"> • Premier Manual Defrost Laboratory Freezer, -30C • 30 Cu. Ft. capacity • One solid swing door, right hinged, not reversible • Five Inner "fast-freeze compartments" with doors for temperature maintenance and recovery • Microprocessor Temperature Controller • Adjustable operating temperature range: -27°-33°C • Standard Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms • Sample simulated digital temperature display; .1° Resolution • Adjustable temperature display in Celsius or Fahrenheit • Audible and visual high/low temperature alarms • Sensor error alarm • Power failure alarm • Resettable 1-point Min/Max History • Remote alarm contacts • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • Four fixed shelves • Keyed door lock • Magnetic door gasket for positive seal • Stainless steel door handle • Caster platform (adds 4" to height) • One probe access port (3/8") • White textured cabinet and door • High density urethane foam cabinet and door insulation

SPEC SHEET ABT-MFP-3030

Product Name:	ABT-MFP-3030
Features:	<ul style="list-style-type: none"> • CFC-free refrigerant (R404a), foam insulation and packaging • Two year parts and labor warranty; plus additional three year compressor parts warranty • Exterior dimensions: 35 ¾" W x 34 ¼" D x 76 ½" H with caster platform; 72 ½" H without caster platform • Manual defrost • 115V, 60 Hz, 10.7 Amps, 1HP • UL/C-UL listed • Shipping weight: 462 lbs.