

SPEC SHEET ABT-MFS-30

Product Name:	ABT-MFS-30
Features:	<ul style="list-style-type: none"> • Standard Manual Defrost Laboratory Freezer • 30 Cu. Ft. Capacity • One solid swing door, right hinged, not reversible • Mechanical thermostat • Adjustable operating temperature range: -15°-25°C • Four fixed shelves • Keyed door lock • Magnetic door gasket for positive seal • Stainless steel door handle • Leveling legs • Probe access port (3/8") • White textured cabinet and door • High density urethane foam cabinet and door insulation • CFC-free refrigerant (R404a), foam insulation and packaging • One year parts and labor warranty; plus additional four year compressor parts warranty • Exterior dimensions: 35 ¾" W x 34 1/4" D x 72 ½" H • Manual defrost • 115V, 60Hz, 6 Amps, 1/3 HP • UL/C-UL listed • Shipping weight: 365 lbs.