

SPEC SHEET CRT-ABT-HC-S12S

| | |
|----------------------|--|
| Product Name: | CRT-ABT-HC-S12S |
| Features: | <ul style="list-style-type: none"> • Upright Controlled Room Temperature Cabinet, Solid Door • 12 Cu. Ft. Capacity • One swing solid door, self-closing (energy efficient, double pane vacuum insulated), right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • Adjustable operating temperature range: 68°F–77°F (20°C – 25°C) • Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable • 2 probes; 1 in air and 1 in sample bottle • Three adjustable shelves • Vaccine storage power cord warning label • Forced draft air circulation • Keyed door lock • Magnetic door gasket for positive seal • Heavy duty door handle • Swivel Casters (3") • LED interior lights are safety shielded and switch controlled • Probe access port (3/8") • Powder coated steel white interior and exterior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) • One year parts and labor warranty, plus an additional four year compressor parts warranty |

SPEC SHEET CRT-ABT-HC-S12S

| | |
|----------------------|---|
| Product Name: | CRT-ABT-HC-S12S |
| Features: | <ul style="list-style-type: none">• Exterior dimensions: 25" W x 29 ¾" D x 65 ¾" H with casters; 62 ¼" H without casters• Cycle defrost• 115V, 60 Hz, 3 Amps, 1/5 HP• ETL listed• Shipping weight: 227 lbs. |