

**SPEC SHEET CRT-ABT-HC-S36S**

<b>Product Name:</b>	CRT-ABT-HC-S36S
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Upright Controlled Room Temperature Cabinet, Solid Door</li> <li>• 36 Cu. Ft. capacity</li> <li>• Two swing solid doors, self-closing, right hinged, not reversible</li> <li>• Microprocessor temperature controller</li> <li>• Digital temperature display</li> <li>• Audible and visual high and low temperature alarms</li> <li>• Remote alarm contacts</li> <li>• Adjustable operating temperature range: 68°F–77°F (20°C – 25°C)</li> <li>• Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual &amp; audible temp alarm</li> <li>• Keyed door locks</li> <li>• Magnetic door gasket for positive seal</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Eight adjustable shelves</li> <li>• Vaccine storage power cord warning label</li> <li>• Forced draft air circulation</li> <li>• Heavy duty door handles</li> <li>• LED interior lights are safety shielded and switch controlled</li> <li>• Swivel Casters (3”)</li> <li>• Probe access port (3/8”)</li> <li>• White powder-coated steel exterior and interior</li> <li>• Hydrocarbon, natural refrigerant (R290)</li> </ul>

**SPEC SHEET CRT-ABT-HC-S36S**

<p><b>Product Name:</b></p>	<p>CRT-ABT-HC-S36S</p>
<p><b>Features:</b></p>	<ul style="list-style-type: none"> <li>• One year parts and labor warranty, plus an additional four year compressor parts warranty</li> <li>• Exterior dimensions: 39 5/8" W x 34 3/4" D x 81 3/4" H with casters; 78 1/4" H without casters</li> <li>• Cycle defrost</li> <li>• 115V, 60 Hz, 4.5 Amps, 1/4 HP</li> <li>• ETL listed</li> <li>• Shipping weight: 360 lbs.</li> </ul>