



- Available w/GL45 screw cap (and drip rings)
- Borosilicate 3.3 hybex glass
- Autoclavable, reusable & breakage resistant
- Permanent white enamel marking area
- Ideal for shaking, media prep or storage

hybex™ Erlenmeyer Flasks are manufactured from high quality, **borosilicate 3.3 glass**. A low coefficient of thermal expansion (3.3), provides temperature resistance that is far superior to that of standard lab glass (or soda lime glass). Hybex glass is designed to withstand high temperatures, making them ideal for shaking applications, media preparation and long-term storage.

The standard Erlenmeyer flasks feature a reinforced rim along with permanent, white enamel graduations for volume measurement and a white identification box for marking or coding.

For storage of media and cultures, the capped flasks include an autoclavable polypropylene cap & drip-ring. Optional cap colors are available, as well as PBT material caps for temperature applications up to 180°C. (Threads are standard GL-45 for compatibility with a wide range of cap types from various manufacturers, including phenolic caps, vented caps, etc.)

Technical Data

Material:	hybex™ borosilicate 3.3 glass
Cap Thread Size:	GL45*
Sterilization:	Autoclave or dry heat
Cap Material:	Polypropylene Autoclavable (up to 140°C)
Included Cap:	Blue

hybex™ Erlenmeyer Flasks

BOROSILICATE 3.3 GLASS (STANDARD, NO CAPS OR THREADS)

ITEM #	DESCRIPTION	QUANTITY
B3110-50	50 mL Erlenmeyer flask	12/pk
B3110-125	125 mL Erlenmeyer flask	12/pk
B3110-250	250 mL Erlenmeyer flask	12/pk
B3110-500	500 mL Erlenmeyer flask	6/pk
B3110-1000	1000 mL Erlenmeyer flask	6/pk
B3110-2000	2000 mL Erlenmeyer flask	1 flask

BOROSILICATE 3.3 GLASS (WITH GL45 CAPS)

ITEM #	DESCRIPTION	QUANTITY
B3220-250	250 mL Erlenmeyer flask with GL45 Cap	12/pk
B3220-500	500 mL Erlenmeyer flask with GL45 Cap	6/pk
B3220-1000	1000 mL Erlenmeyer flask with GL45 Cap	6/pk
B3220-2000	2000 mL Erlenmeyer flask with GL45 Cap	1 flask
B3000-HTC	High temperature red caps, up to 180°C, PBT (GL45)	10/pk
B3000-CAP*	Replacement cap, standard blue (GL45)	10/pk

*Add (-G) for green caps, (-P) for purple caps, (-Y) for yellow caps, (-W) for white caps, or (-AS) for assorted caps (2 of ea. color)



Available in standard (non-threaded)

