

SPEC SHEET PH-ABT-HC-10PS

Product Name:	PH-ABT-HC-10PS
Features:	<ul style="list-style-type: none"> • Premier Pharmacy/Vaccine Compact Refrigerator • 10.5 Cu. Ft. Capacity • One solid swing door, right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Remote alarm contacts • Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) • Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual & audible temp alarm • Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall • Bottom footed shelf to allow air circulation • Back wall evaporator cover • Vaccine storage power cord warning label • Forced draft circulation • Keyed door lock • Magnetic door gasket for positive seal • Leveling legs • LED interior light • Probe access port (3/8") • Powder coat exterior • High density urethane foam cabinet insulation

SPEC SHEET PH-ABT-HC-10PS

Product Name:	PH-ABT-HC-10PS
Features:	<ul style="list-style-type: none"> • HFC-free refrigerant (R600a), foam insulation and packaging • Pharmacy refrigerator/freezer toolkit and temperature logs • Two year parts and labor warranty, plus an additional three year compressor parts warranty • Pyxis®, Omnicell® and AcuDose RX® compatible • Exterior dimensions: 23 ¾" W x 26 ½" D x 59" H • Note: This unit must have 4" clearance on sides and back for adequate ventilation • UL/C-UL listed • ENERGY STAR CERTIFIED • Cycle defrost • 115V, 60 Hz, 1/3 HP, 2.3 Amps • Energy Consumption 2.42 kWh/day • Shipping weight: 183 lbs.