

**SPEC SHEET PH-ABT-HC-S12S**

<b>Product Name:</b>	PH-ABT-HC-S12S
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Standard Pharmacy/Vaccine Solid Door Refrigerator</li> <li>• 12 Cu. Ft. Capacity</li> <li>• One swing solid door, self-closing (energy efficient, double pane vacuum insulated), right hinged, not reversible</li> <li>• Microprocessor temperature controller</li> <li>• Digital temperature display</li> <li>• Audible and visual high and low temperature alarms</li> <li>• Remote alarm contacts</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C)</li> <li>• Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable</li> <li>• Three adjustable shelves (adjustable in ½" increments)</li> <li>• Vaccine storage power cord warning label</li> <li>• Forced draft air circulation</li> <li>• Keyed door lock</li> <li>• Magnetic door gasket for positive seal</li> <li>• Heavy duty door handle</li> <li>• Swivel Casters (3 ½")</li> <li>• LED interior lights are safety shielded and switch controlled</li> <li>• Probe access port (⅜")</li> <li>• Powder coated steel white interior and exterior</li> <li>• High density urethane foam cabinet insulation</li> <li>• Hydrocarbon, natural refrigerant (R290)</li> <li>• Pharmacy refrigerator/freezer toolkit and temperature logs</li> </ul>

**SPEC SHEET PH-ABT-HC-S12S**

<b>Product Name:</b>	PH-ABT-HC-S12S
<b>Features:</b>	<ul style="list-style-type: none"><li>• One year parts and labor warranty, plus an additional four year compressor parts warranty</li><li>• Exterior dimensions: 25" W x 29 ¾" D x 65 ¾" H with casters; 62 ¼" H without casters</li><li>• Cycle defrost</li><li>• 115V, 60 Hz, 3 Amps, 1/5 HP</li><li>• ETL listed</li><li>• Shipping weight: 227 lbs.</li></ul>