

SPEC SHEET PH-ABT-HC-S23S

Product Name:	PH-ABT-HC-S23S
Features:	<ul style="list-style-type: none"> • Standard Pharmacy/Vaccine Solid Door Refrigerator • 23 Cu. Ft. capacity • One swing solid door, self-closing, right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • 2 probes; 1 in air and 1 in sample bottle • Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) • Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual & audible temp alarm • Keyed door lock • Magnetic door gasket for positive seal • Heavy duty door handle • Four adjustable shelves (adjustable in 1/2" increments) • Vaccine storage power cord warning label • Forced draft air circulation • LED interior lights are safety shielded and switch controlled • Swivel Casters (3 1/2") • Probe access port (3/8") • White powder-coated steel exterior and interior • Hydrocarbon, natural refrigerant (R290) • Pharmacy refrigerator/freezer toolkit and temperature logs

SPEC SHEET PH-ABT-HC-S23S

<p>Product Name:</p>	<p>PH-ABT-HC-S23S</p>
<p>Features:</p>	<ul style="list-style-type: none"> • One year parts and labor warranty, plus an additional four year compressor parts warranty • Exterior dimensions: 26 7/8" W x 34 3/4" D x 81 3/4" H with casters; 78 1/4" H without casters • Cycle defrost • 115V, 60 Hz, 3 Amps, 1/5 HP • ETL listed • Shipping weight: 250 lbs.