

SPEC SHEET PH-ABT-HC-S26S

Product Name:	PH-ABT-HC-S26S
Features:	<ul style="list-style-type: none"> • Standard Pharmacy/Vaccine Solid Door Refrigerator • 26 Cu. Ft. capacity • One swing solid door, self-closing, right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • 2 probes; 1 in air and 1 in sample bottle • Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) • Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable • Four adjustable shelves (adjustable in 1/2" increments) • Vaccine storage power cord warning label • Forced draft air circulation • Keyed door lock • Magnetic door gasket for positive seal • Heavy duty door handle • LED interior lights are safety shielded and switch controlled • Swivel Casters (3 1/2") • Probe access port (3/8") • White powder-coated steel exterior and interior • Hydrocarbon, natural refrigerant (R290) • Pharmacy refrigerator/freezer toolkit and temperature logs • One year parts and labor warranty, plus an additional four year compressor parts warranty

SPEC SHEET PH-ABT-HC-S26S

Product Name:	PH-ABT-HC-S26S
Features:	<ul style="list-style-type: none">• Exterior dimensions: 28 3/8" W x 36 1/2" D x 81 3/4" H with casters; 78 1/4" H without casters• Cycle defrost• 115V, 60 Hz, 3 Amps, 1/5 HP• ETL listed• Shipping weight: 275 lbs.