

SPEC SHEET PH-ABT-HC-S36S

Product Name:	PH-ABT-HC-S36S
Features:	<ul style="list-style-type: none"> • Standard Pharmacy/Vaccine Solid Door Refrigerator • 36 Cu. Ft. capacity • Two swing solid doors, self-closing, right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • 2 probes; 1 in air and 1 in sample bottle • Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) • Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual & audible temp alarm • Keyed door locks • Magnetic door gasket for positive seal • Eight adjustable shelves (adjustable in ½" increments) • Vaccine storage power cord warning label • Forced draft air circulation • Heavy duty door handles • LED interior lights are safety shielded and switch controlled • Swivel Casters (3 ½") • Probe access port (¾") • White powder-coated steel exterior and interior • Hydrocarbon, natural refrigerant (R290) • Pharmacy refrigerator/freezer toolkit and temperature logs

SPEC SHEET PH-ABT-HC-S36S

Product Name:	PH-ABT-HC-S36S
Features:	<ul style="list-style-type: none">• One year parts and labor warranty, plus an additional four year compressor parts warranty• Exterior dimensions: 39 5/8" W x 34 3/4" D x 81 3/4" H with casters; 78 1/4" H without casters• Cycle defrost• 115V, 60 Hz, 4.5 Amps, 1/4 HP• ETL listed• Shipping weight: 360 lbs.