

SPEC SHEET PH-ABT-HC-SSP-23FA

Product Name:	PH-ABT-HC-SSP-23FA
Features:	<ul style="list-style-type: none"> • Premier Series Pharmacy Stainless Steel Freezer • 23 Cu. Ft. capacity • One swing solid door, self-closing, right hinged, not reversible • 3 ea. Adjustable Shelves • Digital display microprocessor temperature controller • Adjustable operating temperature: 5°F – -13°F (-15°C - -25°C) • Audible and visual Hi/Lo temperature alarm • 2 probes; 1 in air and 1 in sample bottle • Probe attached to evaporator for defrost termination sensing • Sensor error alarm • Continuous product (simulated) temperature display • Remote alarm contacts • Forced draft air circulation • Keyed door lock • Magnetic door gaskets for positive seal • Integrated door handle • LED interior lights are safety shielded and switch controlled • Probe access port (3/4") • Stainless steel exterior • White-powder coated aluminum interior with 304 stainless steel floor • Easily cleaned smooth interior • Hydrocarbon natural refrigerant (R290) • Two year parts and labor warranty, plus an additional three year compressor parts warranty • Auto defrost

SPEC SHEET PH-ABT-HC-SSP-23FA

Product Name:	PH-ABT-HC-SSP-23FA
Features:	<ul style="list-style-type: none">• Four swivel casters (5 3/4")• Exterior dimensions: 26 7/8" W x 35 3/4" D x 83 1/4" H with casters; 77 1/2" H without casters• 115V, 60 Hz, 8 Amps, 3/4 HP• ETL listed• Shipping weight: 327 lbs.