

Shakers Incubators Centrifuges Balances Rotators Vortexers DryBaths Rockers and MORE!



BeadBath™ Duo Series Laboratory Baths

- Two separate modes: Use with beads or water
- Quik-Cal™ simple in-lab calibration function
- Available in 2L or 10L volumes
- Temperature Ranges:
WaterMode: Ambient +5 to 80°C
BeadMode: Ambient +5 to 70°C



Advantages of Bead Baths:

- Eliminates constant refills (Burn-out protection)
- Easy (accessory free) securing of tubes and bottles
- Allows angled incubations
- Runs constant (no wait for heat-up)



Two Modes:
Choose BEAD or WATER


Item No.	Description	Price
B2402	BeadBath Duo, 2L with lift off lid	\$847
B2410	BeadBath Duo, 10L with hinged lid	\$1,446



BenchMixer™ MP Vortexer

- Accepts a wide range of microplates, including deep well, PCR, ELISA plates, etc.
- Includes tube accessories for vortexing of 1.5/2.0ml, 15ml, or 50ml tubes
- Converts to a traditional "touch" mode vortexer with included cup head

Included!




Item No.	Description	Price
BV1007	BenchMixer™ MP Vortex Mixer	\$499



SureTemp™ Ovens

- Temperatures up to 200°C
- Large viewing window for monitoring
- Stainless steel chamber
- In-lab calibration feature
- Three sizes: 30L, 50L and 150L



Item No.	Description	Price
H2606-30	Benchmark SureTemp Oven 30L	\$1,295
H2606-50	Benchmark SureTemp Oven 50L	\$1,595
H2606-150	Benchmark SureTemp Oven 150L	\$2,895

Temp Range: +5 to 200°C
Stability: +/- 0.5°C at 100°C and +/- 1.0°C at 200°C
Timer: 1 min to 999 hrs or continuous



C2417-R (refrigerated)

MC-24™ Series Microcentrifuges

- Includes COMBI-Rotor™ for tubes and strips
- High speed, up to 16,800xg
- Available with refrigeration
-20°C to +40°C temp. range

Special Offer:
FREE myFuge Mini Centrifuge

Receive a MF8-G w/ purchase of MC24 or C2417-R

*enter MF8-G-VF to your cart for free item



Item No.	Description	Price
MC24	MC-24 Microcentrifuge w/ COMBI-rotor	\$2,693
C2417-R	MC-24 Refrigerated Microcentrifuge w/ COMBI-rotor	\$5,780



Touch Control High Speed



MC24 (ambient)



BeadBug™ and BeadBlaster™ Homogenizers

- Completely lyse, grind & homogenize a variety of sample types
- Available with **REFRIGERATION** to protect delicate samples
- Process up to 24 samples simultaneously
- Wide variety of grinding media in prefilled 2ml & 5ml tubes

Item No.	Description	Price
D1030*	BeadBug 3 Homogenizer	\$998
D1036*	BeadBug 6 Homogenizer	\$2,736
D2400*	BeadBlaster 24 Homogenizer	\$9,460
D2400-R*	BeadBlaster 24R Refrigerated Homogenizer	\$15,098
D1032-05	Prefilled 2ml tubes, 0.5mm zirconium, 50/pk.	\$184
D1032-15	Prefilled 2ml tubes, 1.5mm zirconium, 50/pk.	\$184
D1032-30	Prefilled 2ml tubes, 3mm zirconium, 50/pk.	\$184
D1033-28	Prefilled 2ml tubes, 2.8mm steel, 50/pk.	\$184
D1034-28	Prefilled 5ml tubes, 2.8mm steel, 50/pk.	\$246

Additional grinding media is available in prefilled tubes and loose bulk. Ask for details.



Special Offer:
FREE Prefilled Tube Kit (50pk.)

*Choose 1 free pack of prefilled tubes w/ purchase. Enter tube item # plus *-VF* for free product.



SmartWasher™ 96 Microplate Washer

- Quick set up & easy programming
- Compatible with U, V, C, and flat bottom 96 well plates
- Flexible programming for full plate, partial plate, or strip washing
- Choose 8-channel or 12-channel wash head
- Individually controlled wash & rinse sources

Item No.	Description	Price
MW9600*	SmartWasher™ Microplate Washer for 96 well plates	\$5,237
MW9600-12CM	12-Channel manifold	\$388
MW9600-8CM	8-Channel manifold	\$353


*manifolds not included (bottles included)

SureTherm™ CO₂ Incubators with IncuVIEW™ LCI

- Precision temperature control with 6-sided heating
- Internal microscope option for **Live Cell Imaging**
- Available with:
 - Digital CO₂ and O₂ control
 - High heat decontamination
 - Split inner door

Special Offer

Item No.	Description	Price
H3565-45	ST-45 CO ₂ Incubator with two shelves	\$4,850
H3565-180	ST-180 CO ₂ Incubator with three shelves	\$4,950
H3565-180HD	ST-180 CO ₂ Incubator w/ high heat decontamination	\$5,880
H3565-180HDO2	ST-180 CO ₂ Incubator w/ high heat decontamination split window door & O ₂ control	\$13,101
H3565-45-LCI	IncuView LCI Integrated microscope shelf (14 x 11.5 in.)	\$4,719
H3565-180-LCI	IncuView LCI Integrated microscope shelf (20 x 17 in.)	\$4,916



IncuView LCI

- Live remote monitoring
- Remote adjustment & imaging
- Capture in HD
- Time lapse and video capture
- Adjustable lighting & focus

NEW! View cells without disturbing their CO₂ environment...
On the front display, or remotely from a device!



Roto-Therm™ Plus Incubated Rotator

- Temperature control up to 60°C
- Includes tube holders for 0.5ml to 50ml
- Easy, in-lab temperature calibration
- Available with variable speed & mix modes
 - Oscillating
 - Pausing
 - Rocking

Included Holders

- 41 x 1.5-2ml (2/pk.)
 - 5 x 5-7ml (2/pk.)
 - 4 x 10-15ml (2/pk.)
 - 6 x 50ml* (1ea.)
- *Also accepts 24 x 0.5ml

Temp Range: Ambient +5° to 60°C
Speed Range: 5 to 70rpm

Item No.	Description	Price
RTP2024	RotoTherm™ Plus w/ included tube holders, variable speed	\$907
R2020-PCR	Optional tube holder for 96 x 0.2 ml (8 x PCR strips)	\$77



Incu-Shaker™ Series Shaking Incubators

- **MAGic™** Clamp accessories for **instant-magnetic clamp exchange**
- Available with **CO₂ control & cooling**
- **NEW!** Adjustable clamp for 125, 250 and 500ml flasks

Special Offer:
FREE MAGic Clamp Platform

Receive IS-PLAT or IS-1010-PLAT w/ purchase of Incu-Shaker Shaking Incubator
Enter platform item # plus "-VF" for free product.



H1000-MR-CMB
Adjustable clamp for 125, 250 & 500ml flasks

Incu-Shaker
Mini

Speed Range: **30 to 300rpm**
Temp Range: **Amb. +5° to 70°C**
Platform Dims: **9.5 x 11.5 in.**
Universal Platform: **IS-PLAT**
Flask Capacity:
2 x 1L (IS-1000)
5 x 500 ml (IS-500)
8 x 250 ml (IS-250)
12 x 125 ml (IS-125)

\$3,635



Item: IS-300

Incu-Shaker
10L & 10LR

Speed Range: **30 to 300rpm**
Temp Range: **(H2012, Cooling Model)**
Amb. +5° to 60°C
Amb. -15° to 60°C
Platform Dims: **17.8 x 17.8 in.**
Universal Platform: **IS-1010-PLAT**
Flask Capacity:
5 x 2L (IS-2000)
9 x 1L (IS-1000)
16 x 500 ml (IS-500)
25 x 250 ml (IS-250)
36 x 125 ml (IS-125)

**\$7,570/
\$10,701**



Item: IS-1010 & IS-1010R (Cooling)

Incu-Shaker
Mini CO₂

Speed Range: **30 to 300rpm**
Temp Range: **Amb. +5° to 60°C**
CO₂ Range: **0 to 20%**
Platform Dims: **9.5 x 11.5 in.**
Universal Platform: **IS-PLAT**
Flask Capacity:
4 x 1L (IS-1000)
5 x 500 ml (IS-500)
8 x 250 ml (IS-250)
12 x 125 ml (IS-125)

\$11,120



Item: H3501



Z287-A Microcentrifuge Plus

- Multipurpose, 11 rotor options
- Micro bundle includes 2 rotors!
- Swing out rotors for tubes & plates up to 50ml
- Unique EZ-Scroll™ touch pad
- Compact & economical

Z287-A Micro Bundle Available!

Item No.	Description	Price
Z287-BND-MB	Z287-A Micro bundle with 18x2.0 ml rotor and 2 x microplate rotor	\$4,800

5-Year Warranty!

HERMLE
Benchmark



18 x 1.5/2ml fixed angle (14,000 rpm / 15,994 xg)



Microplate swing out rotor (3,800 rpm / 3,717 xg)



myBlock™ Series Digital Dry Baths

- Digital temperature control up to 105°C
- Hinged lid improves uniformity & conserves energy
- External temperature probe (optional)
- Available with **Quick-Flip™** universal blocks for microtubes (0.2 to 2.0ml), PCR strips & PCR plates

Item No.	Description	Price
LDMB1	myBlock I Single Chamber Dry Bath	\$555
LDMB2	myBlock II Dual Chamber Dry Bath	\$691
LDMB1WB	myBlock I Single Chamber Dry Bath, with 1x LDUB	\$720
LDMB2WB	myBlock II Dual Chamber Dry Bath, with 2x LDUB	\$1,010
BSH-TP1	Optional external temp. probe, for sample temp. control (Blocks sold separately)	\$239

Block Ordering Information

LDUB	Quick-Flip, 24 x 1.5/2.0ml, 14 x 0.5ml and 6 x PCR Strips / 48 x 0.2ml	\$204
BSW02	48 x 0.2ml (6 x PCR Strips)	\$152
BSW1500	24 x 1.5ml	\$152
BSW1520	24 x 1.5/2.0ml	\$152
BSW13	20 x 13mm	\$152
BSW15	12 x 15ml conical	\$152
BSW50	5 x 50ml conical	\$152
BSWPCR1	1 x 96 well PCR plate (for Single Dry Bath only)	\$342

(Additional blocks available, please visit our website for a complete listing)

QuickFlip™ Reversible Block

BSWCMB



24 x 1.5/2.0ml & 14 x 0.5ml

48 x 0.2ml or 6 x PCR Strips*

*Place 2 blocks together in a dual dry bath to accept a 96 well PCR plate



LDUB (QuickFlip)




BSW02



BSW1500



BSW1520




BSW13



BSW15



BSW50



myBlock I (Single Chamber Dry Bath)

Choose one
FREE

With purchase of item below*

(Free with items: IB6000, H3715, B5200, E7500, H5100-HCT, OS-4000)



Value: \$462

Includes COMBI-rotor for tubes & strips



Value: \$335

Speed Range: 200-3200 rpm



Value: \$441

Your choice of:

MyFuge™ 12

Item: MF12-VF

Benchmixer™ V2

Item: VM-1003-VF

RotoMini™

Item: RFP2020-VF

*Promo Code BP1.25

ISO BLOCK™
Two Independent Block Chambers

- Two isolated chambers with independent temperature control
- Digital control up to 105°C
- Hinged lid improves uniformity
- Quick-Flip universal block accepts 24 x 1.5ml, 48 x 0.2ml, 14 x 0.5ml


Additional blocks available for all common tubes & vials

IB6000	isoBlock™ Digital Dry Bath, two independent chambers	\$945
LDUB	Quick-Flip™ Block, 24 x 1.5ml, 48 x 0.2ml, 14 x 0.5ml	\$203
BSW02	Block, 48 x 0.2ml tubes or 6 PCR strips	\$152

HS-10 Plus
Multi-functional Hotplate Stirrer



- Multi-functional design, simultaneously heat & stir up to 5 vessels
- Large, 10 x 10 inch white, ceramic top plate
- Digital speed and temperature control
- Safety LED indicates hot surface



H3715 Benchmark HS-10 Plus Hotplate Stirrer **\$1,049**

SUREAIR™
PCR WORKSTATION WITH UV AND HEPA FILTRATION

- HEPA filter and UV light reduce cross-contamination during PCR preparation
- Compact, 55 x 40 cm footprint
- Cable port allows easy use of centrifuges, vortexers and other instruments
- Includes pipette holder and shelf for storage of pipettes and tips
- Provides clean air (to ISO 5 standards)

B5200	SureAir™ PCR Workstation	\$2,325
Accessories:		
B5200-FIL	SureAir™ Replacement Filter	\$314
B1450-IC	UV-C Indicator Cards, 25 ct.	\$111
B1450-ISS	UV-C Indicator Stickers, 40 ct	\$99

QüadCount™
Automated Cell Counter
for mammalian cells

- Fast, accurate and reproducible counting & live/dead determinations
- Intuitive software with auto-focus and touch-screen interface
- Cost-effective QuadSlides™
- Disposable slides with four individual chambers



E7500	QuadCount Cell Counter	\$4,361
E7500-S1	QuadSlide, 4 chamber slides	\$175

MULTI-THERM™ Touch
COOL - HEAT - SHAKE

- Large color touch screen
- Multiple block options for plates & tubes up to 50ml
- Broad temperature & speed range
- Advanced programming features




Heated lid prevents condensation & evaporation

Speed range: 200-3,000 rpm (block dependent)
Shaker Orbit: 3mm
Temp. range: Ambient -20° to 105°C

H5100-HCT	MultiTherm Touch, Heating/Cooling Vortexer (blocks sold separately)	\$4,231	Available blocks*		
			H5100-05	24 x 0.5ml	\$378
			H5100-MP	Block, microplate	\$378
			H5100-384	Block, 384 well plate	\$648
H5100-HL	Optional heated lid	\$811	*Additional blocks available, see website for details		

ORBI-SHAKER™ CO₂
Remote Controlled, All Environment Shakers

- Designed to fit & maximize capacity in most CO₂ incubators & to withstand temperatures 0° to 65°C
- Horizontally circular 19mm orbit for aeration & mixing
- Adjust settings via the included remote control module
- Instantly exchange flask clamps with, **MAGic Clamp™** accessories



Speed range: 30-300 rpm
Max capacity: 2 x 2L or 5 x 1L flasks (BT4001)
5 x 2L or 9 x 1L flasks (BT4011)

OS-4000	Orbi-Shaker CO ₂ Shaker with remote control, includes 13x12" non-slip platform	\$3,879
OS-4000-XL	Orbi-Shaker CO ₂ XL Shaker with remote control, includes 6x16.5" non-slip platform	\$5,379

Free offer with purchase of designated items at listed price. Cannot be combined with other discounts or pricing. Free item must be of equal or lesser value than the lowest-priced item. Not valid on prior purchases. Special price applies only to the specified item; no substitutions. Promo valid until 04/30/2025 or while supplies last. All terms and conditions apply.



ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998