

Gobble up these amazing BLACK FRIDAY DEALS!


**BLACK
FRIDAY
SALE**
THREE WEEKS
LONG



myFUGE™ 12

- * COMBI-Rotor™ for tubes and strips
 - 12 x 1.5/2.0 mL tubes
 - 32 x 0.2 mL PCR tubes
 - 4 x PCR strips (8 x 0.2 mL)
- * Twice the capacity of most mini centrifuges
- * Near silent operation
- * Starts and stops by closing and opening the lid

MF12 myFuge™ 12 Mini Centrifuge with COMBI-Rotor \$425.10/2
C1008-A5-6 Adapter pack, 0.5 mL, pack of 6 \$24.08




BUY ONE AT \$425.10 - GET ONE FREE!

Shake things up!

BEADBUG
MICROTUBE HOMOGENIZER

- * Powerful mixing action for lysis, grinding and homogenization
- * Eliminates cross contamination
- * Available in 3 or 6 tube capacity
- * Cold room safe
- * Wide variety of grinding media available



D1030 BeadBug 3, 3 x 2 mL tube capacity \$925.28
D1036 BeadBug 6, 6 x 2 mL tube capacity \$2457.95
MF8* myFuge Mini 6,000rpm/2,000xg, 8 place rotor ~~\$324.21~~ FREE!
*Add (-B) for Blue, (-G) for Green, (-P) for Purple, (-R) for Red, (-C) for clear
D1032-30 TriplePure Zirconium Beads, 3 mm, pk 50 tubes \$166.69
Additional bead materials/sizes available



Call Us Today!
BLACK FRIDAY
DEALS END
12/2/2022

FREE MYFUGE™ MINI WITH PURCHASE - \$324.21 VALUE!

Hot stuff!

Digital Hot Plate Stirrer

- * Large 7 x 7" ceramic surface
- * 5 L stirring capacity
- * Digital speed and temperature control
- * Optional probe for controlling by sample temp

H3770-HS 7 x 7" Digital Hot Plate Stirrer ~~\$436.91~~ \$339.00
H3770-TP Probe for sample temp control \$126.00
H3770-CS Clamp set, includes rod and clamp \$87.36



HUGE DISCOUNT! SAVE \$97.91!



ORDER TODAY!

midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998

Black Friday BLOWOUT Deals



BenchBlotter™ Benchmark Platform Rocker

- * Perfect for blot washing and gel staining
- * Large 12 x 12 inch platform
- * Gentle mixing - 12 rpm, $\pm 8^\circ$ tilt
- * Cold room and incubator safe
- * Compact, low profile design



BR1000
BR1000-DIMPLED
BR1000-STACK

BenchBlotter™ Platform Rocker w/non-slip mat
Optional dimpled mat for tubes
Stacking platform

\$528.65/2
\$ 43.40
\$124.60

BUY ONE AT \$528.65 - GET ONE FREE!



SMARTSEAL™ Semi-Automated Plate Sealer

- * Fast, reproducible sealing of well plates
- * Adjustable time, temperature and pressure
- * Simply load plate and press run to seal
- * Accepts a wide range of plate types and heights:
Deep well storage plates, assay plates, PCR plates

**AMAZING
BLACK
FRIDAY
DEAL!

SAVE
\$1000!**



MS1000
Sealing Films
MS1000-F1
MS1000-F2
MS1000-PCR1
MS1000-PCR2
MS1000-PCR3

SmartSeal™ Semi-Automated Plate Sealer

~~\$4695.00~~ \$3695.00

Foil, peelable for PP plates, 100/pk	\$79.00
Foil, strong bond, for cold storage & DMSO, PP plates, 100/pk	\$139.00
Optically clear, for qPCR, peelable, for most plates, 100/pk	\$139.00
Optically clear, for qPCR, strong bond, for PP plates, 100/pk	\$139.00
Optically clear, for qPCR, pierceable, for most plates, 100/pk	\$139.00

LIMITED TIME OFFER! - SAVE \$1000!



MULTI-THERM™ Touch COOL - HEAT - SHAKE

- * Large touch screen for simple and easy programming and viewing of status
- * Advanced programming features include adjustable ramp rates, preheat and clockwise/ counterclockwise mixing
- * Multiple block options for plates and tubes up to 50 mL
- * Broad temperature and speed range to meet a variety of applications



MTHC-1500	MultiTherm™ Touch Shaker w/heating/cooling	\$3949.00
MTB-CMB	Combination block: 15 x 0.5 mL & 20 x 1.5mL	\$350.00 FREE!
H5100-HL	Heated lid for MultiTherm Touch	\$750.00 FREE!

*Blocks sold separately. Additional blocks available. See www.midsci.com.



FREE BLOCK AND HEATED LID WITH PURCHASE - \$1100 VALUE!

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998

More Black Friday BLOWOUT Deals

SmartReader™ 96 Microplate Absorbance Reader

- * Stand alone operation - no computer needed
- * User-friendly, touchscreen programming
- * Integrated plate shaker
- * Supplied with 405, 450, 492 and 630nm filters
Additional filters sold separately



MR9600-MID	SmartReader™ 96 Microplate Reader	\$6050.00	\$5050.00
MR9600-TP	Thermal Printer for MR9600	\$341.25	
MR9600-PA	Extra paper for printer, pk of 4	\$32.55	


ACCURIS
INSTRUMENTS

BLACK FRIDAY SALE! SAVE \$1000!



DryCHILL™ by Benchmark

- * Ice-free, stable cooling for sample tubes
- * Unique design reduces condensation
- * Maintains 4°C for up to +200 minutes
- * 4 sizes for various tubes up to 50 mL and PCR plates
- * Visual temperature indicator



DC9602	DryCHILL™ Cooler, 96 x 0.2 mL tubes/PCR plate	\$395.00
DC4020	DryCHILL™ Cooler, 40 x 1.5/2 mL tubes	\$395.00
DC0650	DryCHILL™ Cooler, 6 x 50 mL tubes	\$445.00
DC1215	DryCHILL™ Cooler, 12 x 15 mL tubes	\$445.00
VM1*	Vornado Mini Vortexer, gray cup*	\$211.89 FREE

*Other cup colors available - use -B for blue, -G for green, -P for purple and -R for red

FREE VORNADO™ MINI VORTEXER! \$211 VALUE!

PLATEFUGE™ MicroPlate MicroCentrifuge

- * Quickly spins down droplets in PCR plates
- * Unique, patented swing-out rotor design prevents spillage
- * Half the size of traditional centrifuges
- * Accepts all popular types of PCR plates (skirted, semi, and non) and most microtiter plates



LOADING RUNNING

SAVE \$217!

PCRCENT	PlateFuge™ Microcentrifuge	\$812.05	\$595.00
PCRCENT-A02	Adapter for 0.2ml tubes and strips	\$59.95	
PCRCENT-DWMP	Optional rotor for deep well plates	\$206.01	

ASPIRE LABORATORY ASPIRATOR

- * Compact design with space saving base
- * Built-in vacuum pump with adjustable pressure and vacuum gauge
- * Safe and precise aspiration
- * Supplied with 2 L PC collection bottle, antibacterial filter, tubing, handheld vacuum controller, adapter for 200 µL tips and 40 mm fine tip steel adapter



BIG SALE

SAVE \$300!



V0020	Aspire Laboratory Aspirator	\$968.00	\$668.00
-------	-----------------------------	---------------------	----------

Additional accessories available, including multichannel adapters.

Check out these amazing prices!



myBLOCK™ HL
With Heated Lid

CHOOSE A FREE BLOCK!



- * Heated lid prevents evaporation and condensation
- * Compact design - fits in the palm of your hand
- * Digital control of block and lid temperature
- * Fast heating - reaches 37°C in 2 minutes

Perfect for enzymatic reactions, LAMP assays and other isothermal reactions

MB100-HL	myBlock™ HL Dry Bath with Heated Lid*	\$465.38	
	Choose ONE of the blocks below, FREE		
MB100-02	Block, 40 x 0.2 mL or 5 PCR strips	\$100.24	FREE
MB100-05	Block, 24 x 0.5 mL	\$100.24	FREE
MB100-15HL	Block, 15 x 1.5 mL	\$100.24	FREE

*Blocks sold separately, offer includes one free block

FREE BLOCK WITH PURCHASE - \$100 VALUE!

BenchMixer™ V2



- * Reversible top: use "cup" or "flat/dimpled" surface
 - * Powerful motor for instant vortexing
 - * Q-Drive™ dynamic balancing system minimizes vibration
 - * Variable speed with continuous or touch operation
- BOGO - Put one on every bench!**

VM1003 BenchMixer V2 with reversible top \$307.26/2



BUY ONE AT \$307.26 - GET ONE FREE!

This sale carries weight!


Accuris™ Tx Semi-Micro Balance

- * Upgraded weighing performance with
- * Electromagnetic Force Restoration (EFR) technology
- * 120 g capacity
- * Semi-micro class - readability of 0.00001 g up to 62 g, and 0.0001 g up to 120 g
- * Color touch screen user interface
- * RS232 for data output



W3002A-120 Analytical Balance Tx, internal calibration, 120g ~~\$3466.20~~ \$2466.00

SAVE \$1000 ON A SEMI-MICRO BALANCE!



ACCURIS
INSTRUMENTS

Free offers require purchase of designated items at listed price. Offers expire 12/02/22. Not valid on previous purchases. Special price applies only to item specified. Other capacities or items cannot be substituted. Use promo code TURKEY22.

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998