



Gobble up these amazing BLACK FRIDAY DEALS!

BLACK FRIDAY SALE

ALL MONTH
LONG

BENCHMIXER™ MP

Vortex Mixer for Tubes and Plates

NEW



Supplied
with unit

- * Exceptionally versatile, supplied with:
 - Microplate head (fits a variety of plates)
 - Adapters for 1.5/2.0mL, 15mL and 50mL tubes
 - Cup head to convert to standard vortexer
- * Digital control
 - 200-3,000rpm, up to 999 minutes or continuous

BV1007 BenchMixer MP for microplates and tubes with adapters for 30x1.5/2mL, 12x15mL, 4x50mL tubes and cup head ~~\$499~~ \$349

INTRODUCTORY PRICE - SAVE \$150!



MC-24R™

REFRIGERATED MICROCENTRIFUGE

- * Touch control and digital display
- * Powerful refrigeration, as low as -20°C
- * High speed - 16,800 xg
- * Quiet operation
- * Supplied complete with COMBI-Rotor™
 - Accepts 24x1.5/2.0mL tubes AND 2xPCR strips



COMBI-Rotor

C2417-R MC-24R Refrigerated Microcentrifuge ~~\$5,695.99~~ \$4,195

Optional rotors:

C2417-2420H 24x1.5/2.0mL hermetically sealed rotor \$841.40

C2417-3202 4xPCR strip rotor (32x0.2mL) \$795.00

SPECIAL BLACK FRIDAY PRICE - SAVE \$1500!

Call Us Today!
BLACK FRIDAY
DEALS END
12/6/2024



TC 9639 Gradient Thermal Cycler

- * Intuitive touch screen control
- * Program wizard - quick, easy and foolproof
- * Multifformat block for 96x0.2mL tubes, PCR strips, 39x0.5mL tubes or a 96 well plate
- * Advanced programming features such as adjustable ramp rate and time/temperature increments/decrements



SALE!
FREE
PlateFuge
with
Purchase!



BSTC-96 TC 9639 Gradient Thermal Cycler ~~\$5995~~ \$4,995
PCRCENT-VF PlateFuge™ Microplate Microcentrifuge ~~\$894.54~~ **FREE!**

FREE PLATEFUGE - \$894 VALUE PLUS \$1000 DISCOUNT!



ORDER TODAY!

midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998

MIDSCI™
Laboratory Equipment & Supplies



Black Friday BOGO Deals




SMART SAFE Nucleic Acid Stain GLOW

- * Excellent replacement for Ethidium Bromide
- * Excitation by UV or blue light
- * Non-mutagenic, environmentally safe
- * Available as a prestain (PS) or loading dye (LD)

Buy one, get one FREE! Mix or Match!
Stock up and SAVE!

BE4500-PS SmartGlow Prestain, 1mL (20,000X) \$124.95
BE4500-LD SmartGlow Loading Dye with stain, 1mL (6X) \$124.95

For free product: Use item number plus "-VF"



BUY ONE AT \$124.95 - GET ONE FREE!



BEADBUG MICROTUBE HOMOGENIZER

- * Powerful mixing action for lysis, grinding and homogenization
- * Eliminates cross contamination
- * Personal sized - 3 tube capacity
- * Cold room compatible
- * Wide variety of grinding media available

D1030 BeadBug 3, 3x2mL tube capacity \$998.00
D1030-VF PROMO ITEM! BeadBug 3 ~~\$998.00~~ FREE!
D1032-30 TriplePure Zirconium Beads, 3mm, pk 50 tubes \$182.86

Additional bead materials/sizes available - see special offer below!

Perfect for small labs and student labs!



BUY ONE AT \$998 - GET ONE FREE!



Benchmark Homogenizer Bead Tubes

- * Convenient, prefilled, 2mL bead tubes
- * Compatible with most brands of bead beaters
- * Variety of bead sizes and materials to meet a range of sample needs
- * TriplePure™ line certified nuclease, protease and DNA free

| | | |
|------------------|--|-----------------|
| D1032-30 | TriplePure Bead Tubes, Prefilled, 3.0mm zirconium, pk 50 tubes | \$182.86 |
| D1032-15 | TriplePure Bead Tubes, Prefilled, 1.5mm zirconium, pk 50 tubes | \$182.86 |
| D1032-10 | TriplePure Bead Tubes, Prefilled, 1.0mm zirconium, pk 50 tubes | \$182.86 |
| D1032-05 | TriplePure Bead Tubes, Prefilled, 0.5mm zirconium, pk 50 tubes | \$182.86 |
| D1032-01 | TriplePure Bead Tubes, Prefilled, 0.1mm zirconium, pk 50 tubes | \$182.86 |
| D1032-SK | TriplePure Starter Kit, 10 each of 0.1, 0.5, 1.0, 1.5 and 3.0mm | \$182.86 |
| D1031-10 | Prefilled Bead Tubes, 1.0mm glass, pk of 50 tubes | \$131.00 |
| D1031-05 | Prefilled Bead Tubes, 0.5mm glass, pk of 50 tubes | \$131.00 |
| D1031-01 | Prefilled Bead Tubes, 0.1mm glass, pk of 50 tubes | \$131.00 |
| D1033-28 | Prefilled Bead Tubes, 2.8mm steel, pk of 50 tubes | \$182.86 |
| D1033-30G | Prefilled Bead Tubes, Garnet Shards, 1x6mm zirconium, pk/50 | \$182.86 |

For free product: Use item number plus "-VF"



BUY 2 PACKS, GET ONE FREE - MIX OR MATCH!

BOGO

More BOGOs and Deals



EZ Pack Agarose Tablets

- * 0.5g tablets eliminate weighing
- * Fast dissolving - just 2 minutes
- * Popular Benchmark Agarose LE
- * Convenient blister pack for safe, easy dispensing

A2501 EZ Pack Agarose, 200 tablets (100g) \$194.16
 A2501-VF Promo Item! EZ Pack Agarose ~~\$194.16~~ FREE!



BUY ONE FOR \$194, GET ONE FREE!


Spin microtubes and PCR strips without changing the rotor!



myFUGE™ 12

- * COMBI-Rotor™ for tubes and strips
 - 12x1.5/2.0mL tubes
 - 32x0.2mL PCR tubes
 - 4xPCR strips (8x0.2mL)
- * twice the capacity of most mini centrifuges
- * Near silent operation
- * Starts and stops with closing and opening the lid

MF12 myFuge 12 Mini Centrifuge with COMBI-rotor \$457.47
 MF12-VF Promo Item! myFuge 12 Mini Centrifuge ~~\$457.47~~ FREE!
 C1008-A5-6 Adapters for 0.5mL tubes, pk of 6 \$26.99



COMBI-Rotor

BUY ONE, GET ONE FREE!

BUY ONE AT \$457, GET ONE FREE!

BactiZapper™

- * Replaces open flame burners for sterilization of loops, needles and tubes
- * Adjustable tilt for ergonomic placement
- * No splatter
- * Optional slide dryer and loop holder



BUY 1 GET 1 FREE
GET 1 FREE
BUY 1 GET 1 FREE

BZT-1001 BactiZapper Tilt Micro-Sterilizer \$576.59
 BZT-1001-VF Promo Item! BactiZapper ~~\$576.59~~ FREE!
 BZ-1001 Loop holder attachment \$71.38
 BZ-1003 Slider dryer attachment for 3 slides \$141.42

BUY ONE BACTIZAPPER - GET ONE FREE!

Roto-Mini™ Plus

- * Variable speed rotation, 5-70rpm
- * Advanced shaking/mixing modes: Mix (Oscillating), Pause, Rock
- * Includes holders for common tubes:
 - 11 x 1.5/0.2mL (pack of 2)
 - 5 x 10-15mL (pack of 2)
 - 5 x 5-7mL (pack of 2)
 - 6 x 50mL (pack of 2)



RMP2024 Roto-Mini Plus Variable Speed Rotator \$530.71
 RMP2024-VF Promo item! Roto-Mini Plus ~~\$530.71~~ FREE!
 R2020-RP Mini rocking platform w/rubber mat, 4.5x3 in. \$114.43
 R2020-PCR Holder, 96x0.2mL, 12xPCR strips, 1xPCR plate \$88.07

BUY ONE, GET ONE FREE - \$530 VALUE!

Check out these UNBELIEVABLE DEALS!

ACCURIS Compact Balances



- * Ideal for basic weighing needs
- * Multiple modes: weigh, count, percent
- * LCD display, 0.01g readability
- * Gold-plated, ceramic capacitors
- * Easy, in-lab calibration

Buy one, get one FREE! Mix or match sizes!*

| | | |
|------------------|--|-----------------|
| W3300-120 | Accuris Compact Balance, 120g capacity | \$327.22 |
| W3300-300 | Accuris Compact Balance, 300g capacity | \$381.68 |
| W3300-500 | Accuris Compact Balance, 500g capacity | \$381.68 |

**Least expensive item is free, use item number plus "-VF" for free product*

BOGO - MIX or MATCH - SAVE UP TO \$381

**BOGO
FREE!**



**FREE
myFUGE
Mini!**



UV CLAVE Ultraviolet Chamber

- * Improves quality of PCR research by reducing potential contaminants
- * Multiple lamps and reflective chamber walls provide true 360° irradiation
- * Proprietary UVC transparent shelf, eliminates potential for missed coverage (blind spots)
- * Inactivates DNA, RNA, viruses & bacteria on PCR vessels
- * Adjustable cycle time



| | | |
|-----------------|--|-----------------------|
| B1450 | UVClave UltraViolet Chamber (research use only) | \$573.92 |
| MF8-B-VF | myFuge Mini Centrifuge with 1.5/2mL tube rotor, 0.2mL strip tube rotor and 0.5mL tube adapters | \$346.90 FREE! |

FREE myFUGE MINI WITH PURCHASE - \$346 VALUE!

PLATEFUGE™ MICROPLATE MICROCENTRIFUGE

- * Quickly spins down droplets in PCR plates
- * Unique, patented swing-out rotor design prevents spillage
- * Half the size of traditional centrifuges
- * Accepts all popular types of PCR plates (skirted, semi, and non) and most microtiter plates



| | | |
|---------------------|--------------------------------------|----------------------------------|
| PCRCENT | PlateFuge Microplate Microcentrifuge | \$894.54 |
| PCRCENT-VF | Promo Item! PlateFuge | \$894.54 FREE! |
| PCRCENT-A02 | Adapter for 0.2mL tubes and strips | \$67.26 |
| PCRCENT-DWMP | Optional rotor for deepwell plates | \$231.34 |



BUY ONE, GET ONE - \$894 VALUE!

Free offers require purchase of designated items at listed price. Offers are valid 11/04/2024 to 12/06/2024. Not valid on previous purchases. Special price applies only to item specified. Other capacities or items cannot be substituted. Use promo code TURKEY24

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998

MIDSCI™
Laboratory Equipment & Supplies