

EZcap™

1x8 PCR

FilmStrips™

For sealing a PCR 8-tube strip or a single row of a 96-well plate

i Application:
PCR and short-term storage

🌡 Temperature Range:
-20°C to +120°C



| CATALOG # | DESCRIPTION | QUANTITY |
|--------------|--------------------------|------------|
| US-CS1X8-400 | EZcap 1x8 PCR FilmStrips | 400 Strips |

Features and Benefits

- Certified DNase- & RNase- free
- Features an extra strong, thick adhesive to minimize evaporation
- Each sheet contains 8 strips on a continuous liner
- An economical and more user-friendly alternative to PCR strip-caps
- Precision cut to a smaller dimension for sealing one PCR 8-tube strip or a single row of a 96-well plate whenever rows must be selectively protected or accessed
- Remove easily with no splash back as can occur with strip-caps

Bundling Opportunities

- PCR 8-tube strips
- PCR plates
- Deep well/ storage plates
- PCR reagents and buffers
- Thermal cyclers
- Pipettors and pipet tips
- Reservoirs

