

NEW


Economical Vertical Chest Freezer Racks for 2" Boxes

Engineered for Safety.
Optimized for Savings.

- Lightweight alternative to conventional Vertical Stainless Steel Freezer Racks
- Constructed of corrosive resistant stainless steel
- Double Welds, Smooth Edges
- Handle for easy rack retrieval from freezer/tank
- Locking Rod included to hold in boxes
- Holds standard 2" H Cryo Boxes
- Lifetime Warranty

Economical Vertical Freezer Racks protect and organize samples in Chest Freezers and Liquid Nitrogen Tanks.

Easy, safe and convenient sample retrieval with attached handle, with locking rod included to hold boxes in the rack. These durable racks offer a lifetime warranty.



Specifications

	Description	Capacity	D x W x H
CFE-10-2	Economical Vertical Freezer Rack for 2" Boxes, Holds 10 Boxes	10 - 2 in. Boxes	140 x 143 x 566 mm 5-1/2 x 5-5/8 x 22-5/16 in.
CFE-11-2	Economical Vertical Freezer Rack for 2" Boxes, Holds 11 Boxes	11 - 2 in. Boxes	140 x 143 x 622 mm 5-1/2 x 5-5/8 x 24-1/2 in.
CFE-12-2	Economical Vertical Freezer Rack for 2" Boxes, Holds 12 Boxes	12 - 2 in. Boxes	140 x 143 x 678 mm 5-1/2 x 5-5/8 x 26-5/8 in.
CFE-13-2	Economical Vertical Freezer Rack for 2" Boxes, Holds 13 Boxes	13 - 2 in. Boxes	140 x 143 x 724 mm 5-1/2 x 5-5/8 x 28-1/2 in.
CFE-15-2	Economical Vertical Freezer Rack for 2" Boxes, Holds 15 Boxes	15 - 2 in. Boxes	140 x 143 x 845 mm 5-1/2 x 5-5/8 x 33-1/4 in.

SEPT 2025