Pharmacy Refrigerator HYC-118





LED Display

Microprocessor control, digital display and temperature adjustment with an increment of 0.1℃



Multiple Alarm Functions

Capable of alerting failures due to high and low temperature, sensor error, power failure, low battery, door ajar and remote alarm



Ergonomic Design

810mm cabinet height makes it easy to place under test bench, saving storage space.



Lockable Design

Safety lock to prevent unauthorized access



Better Cabinet Temperature Uniformity, Auto Defrost

The test point of performance in the cabinet increased from 2 points to 9 points.



Energy Efficient Hydrocarbon (HC) Refrigeration System

Optimised refrigeration system, utilising high-efficiency HC compressor significantly reduces energy consumption.



Greener by Design

HC refrigerants are natural non-toxic and more environmentally friendly.



Specifications ()

HYC-118

| Technical Data | |
|----------------|---------------------|
| Cabinet Type | Upright, solid door |
| Climate Class | N |
| Cooling Type | Forced air cooling |
| Defrost Mode | Auto |
| Refrigerant | HC |
| Noise ((dB(A)) | 41 |

| Performance | |
|-----------------|-----|
| Temp Range (°C) | 2~8 |



| Control | |
|------------|----------------|
| Controller | Microprocessor |
| Display | LED |

| Electrical Data | | |
|-----------------------|------------|--------|
| Power Supply(V/HZ) | 220-240/50 | 115/60 |
| Power(w) | 190 | 190 |
| Electrical Current(A) | 1.4 | 1.4 |

| Functions | |
|---------------|---|
| High/Low Temp | Y |
| Remote System | Υ |
| Power Failure | Y |
| Sensor Error | Υ |
| Low Battery | Y |
| Door Ajar | Υ |
| Remote System | Υ |
| | |

| Accessories | |
|-----------------|-----|
| Caster | Υ |
| Foot | Υ |
| Porthole | Υ |
| Shelves/Drawers | 3/1 |

| Dimensions | | |
|------------------------------|---------|----------------|
| Capacity(L/Cu.Ft) | 118/4.2 | |
| Net/Gross | (kg) | 41/46 |
| Weight(approx) | (lbs) | 90.5/101.5 |
| Interior Dimension(w*D*H) | (mm) | 515*415*630 |
| | (in) | 20.3*16.3*24.8 |
| Exterior | (mm) | 597*635*810 |
| Dimension(W*D*H) | (in) | 23.5*25.0*31.9 |
| Packing Dimension(W*D*H) | (mm) | 680*690*910 |
| | (in) | 26.8*27.2*35.8 |
| | | |

| Other | | |
|-------------|----|----|
| Certificate | CE | UL |