## Intelligent Protection of Life Science

# Pharmacy Refrigerator **HYC-290**





### **Frequency Conversion Energy Saving**

Variable frequency compressor with optimized cooling system, more energy efficient.



#### **Environmentally Friendly**

Hydrocarbon compressor and refrigerant, the refrigerant does not contain any fluorine or chlorine which is super friendly to the environment.



#### Large Space, Superior Uniformity

The space is efficiently used, and the utilization rate is high and the temperature uniformity is ±3°C.



#### **Low Noise**

Optimized system and noise reduction design to reduce the noise level to 41dB(A) (220V/50Hz).



#### Stable performance, multiple alarm modes

Three alarm modes ( light ,sound ,remote); six types of alarm



**Microprocessor Control System** 



Self Closing Door with 90° Stay Open Feature.



**Precise Temperature Control** 



**Security lock** 



**Optional USB and Printer** 



# **Scope of Application**



This product is suitable for pharmacy, pharmaceutical manufacturers, anti-epidemic stations, health centers and hospitals as key examples, which is used to store biological products and articles that need to be stored at  $2\sim8$ °C.

# **Specifications**



## **HYC-290**

Technical Data		
Cabinet Type	Upright, glass door	
Climate Class	N	
Cooling Type	Forced air cooling	
Defrost Mode	Auto	
Refrigerant	HC	
Noise ((dB(A))	43	43

Performance	
Temp Range(°C)	2~8

Control	
Controller	Microprocessor
Display	LED

Electrical Data		
Power Supply (V/Hz)	220~240/50&60	115/60
Power (W)	350	200
Electrical Current (A)	1.8	2.2

Dimensions		
Capacity (L/Cu.Ft)	290/10.2	
Net/Gross	(kg)	105/127
Weight (approx)	(lbs)	231.6/280.1
Interior Dimension (W*D*H)	(mm)	530*555*1080
	(in)	20.9*21.9*42.5
Exterior Dimension (W*D*H)	(mm)	665*710*1665
	(in)	26.2*28.0*65.6
Packing Dimension (w*D*H)	(mm)	690*790*1810
	(in)	27.1*31.0*71.1
Container Load (20'/40'/40'H)	21/42/42	



Accessories		
Caster	Y	
Foot	Y	
Porthole	Y	
Shelves/Drawers	5/-	-
USB Interface	Optional	Υ
Temperature Recorder	Optional	

Alarms	
High/Low Temp	Υ
Remote System	Y
Power Failure	Υ
Sensor Error	Υ
Low Battery	Υ
Door Ajar	Υ

Other		
Certificate	CE	UL