

Free Laser Engraving HALO

MidSci™
Laboratory Equipment & Supplies

ERGONOMIC DESIGN

Lightweight with stress point reducing contoured edges and soft-ejection, Halo pipettes are designed for comfort.

ULTRA LIGHT WEIGHT

Up to 40% lighter than the market leaders, HALO™ pipettors are among the lightest ever produced.

SIMPLE CALIBRATION

In-lab calibration is made possible with the calibration key supplied with every pipette.

SUPERIOR ACCURACY

With a 4 digit micrometer, and overlapping volume ranges, aspirations of even fractional volumes are achievable with repeatable precision.

SUPERIOR ACCURACY

Feel confident knowing that your pipette is built to last and covered by a three year manufacturer's warranty.




HALO™ Starter Kit

For a limited time, order a HALO™ Starter Kit #H6800-SK2 and receive free laser-etched personalization on each pipettor. Choose up to 40 custom characters to be permanently engraved on your pipettes.*

Item #	Description	Price
H6800-SK2	HALO Starter Kit, Includes H6800-10, H6800-20, H6800-200, H6800-1M and Carousel Stand, with FREE* Personalization on Pipettes	\$699.00

Included in Kit

Item #	Volume Range	Accuracy	Reproducibility	Price
H6800-10	1 – 10µL	±2.5 – 1.0%	≤1.50 – 0.80%	\$249.00
H6800-20	2 – 20µL	±3.0 – 0.9%	≤2.00 – 0.40%	\$249.00
H6800-200	20 – 200µL	±3.0 – 0.6%	≤1.00 – 0.15%	\$249.00
H6800-1M	100 – 1000µL	±3.0 – 0.6%	≤0.60 – 0.20%	\$249.00
H6800-CR	6-Place Carousel Stand			\$89.50



*Starter kit includes a 10µL, 20µL, 200µL, and 1000µL pipettors and a 6-place carousel stand. 40 custom characters to be spread out over 2 lines of 20. Offer Ends 12/31/23.

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998