

## Centrifuge Tubes: 5mL to 500mL

- Vivid graduations
- A perfect size for any application
- Marking area on most sizes
- RNase, DNase, & Pyrogen free
- 25mL tube fits 50mL rotor

Buy 3 Centrifuge tubes and receive a  
**FREE TUBE RACK**

Buy 3 Centrifuge tubes and receive a  
**FREE TUBE RACK**



5mL  
C1005-SC5-R

10mL  
C2410

15mL  
C2600

25mL  
C2625

30mL  
MID-2630

50mL  
C2602


50mL  
C2605

250mL  
C1250

500mL  
C1500

Specifications		
Max RCF*:	5mL	25,000xg
	10mL	12,000xg
	15mL	12,000xg
	25mL	17,000xg
	30mL	7,000xg
	50mL Conical	15,000xg
	50mL Self-Standing	7,000xg
	250mL	6,000xg
	500mL	6,000xg
Temp Range:	All tubes: 86°C to 121°C	
Material:	Tubes: Polypropylene	
	Caps: Polyethylene	

\* refers to liquids w/ a specific gravity of 1.0, using appropriate rotors, cushions, and adapters.



C2592

R1070

Item	Description	Qty	Price
C1005-SC5-B	Non-sterile 5mL MacroTubes® w/ screw caps packed separately	500	\$99.00
C2535	Sterile 5mL MacroTubes® w/ screw caps, packed in bags of 100	500	\$105.00
C1005-SC5-R	Sterile 5mL MacroTubes® w/ attached screw caps in foam racks	500	\$112.35
C2410	10mL PP (16x104mm), flat screw cap, 50 tubes per sterile bag	1000	\$234.00
C2600	15mL PP (17x118mm), flat screw cap, 50/foam rack in sterile bag	500	\$172.50
C2601	15mL PP (17x118mm), flat screw cap, 50/sterile bag	500	\$162.00
C2592-VF	Rack for 50 x 5mL, 10mL, 15mL tubes, flip-top or screw top	1	<del>\$47.40</del> <b>FREE*</b>

\* With purchase of 3 centrifuge tubes

Buy 3 packs of tubes and receive a  
**FREE TUBE RACK**

C2592-VF (50 x 5mL) or R1070-VF (30 x 15mL & 20 x 50mL)

Item	Description	Qty	Price
C2625-B	25mL PP (29x77mm), with screw cap, sterile, 8 bags of 25 tubes	200	\$101.20
C2625-R	25mL PP (29x77mm), with screw cap, sterile, 8 foam racks of 25 tubes	200	\$104.65
MID-2630	30mL PP, free standing, flat screw cap, 25/sterile bag	500	\$169.05
C2602	50mL PP (29x115mm), flat screw cap, 25/foam rack in sterile bag	500	\$199.50
C2603	50mL PP (29x115mm), flat screw cap, 25/sterile bag	500	\$179.40
C2605	50mL PP, free standing, flat screw cap, 25/sterile bag	500	\$205.85
R1070-VF	Rack for 30x15 & 20x 25mL or 50mL (17x20x6cm), Orange	5	<del>\$69.19</del> <b>FREE*</b>

\* With purchase of 3 centrifuge tubes

Item	Description	Qty	Price
C1250	250mL PP (59.8x160mm), with screw-cap, 14 bags of 5 tubes	70	\$234.00
C1500	500mL conical centrifuge tube, PP, with PE plug type screw cap. Packed 5 per sterile bag, 40 tubes (8 bags) per case.	40	\$239.00

## Cryogenic Vials

- Temperatures as low as -196°C
- Self-standing star foot
- Leakproof silicone O-ring
- Resealable bags of 50
- Gamma sterilized
- Free locking rack
- RNase, DNase, pyrogen, & endotoxin free




**FREE Locking Rack**  
in Every Case

Choose External Thread (first column) or Internal Thread (second column):

External	Internal	Description	Qty	Price
V3801	V4801	Cryogenic Vials, With silicone o-ring, 1mL, sterile, printed graduations and marking area, self-standing star foot, 50-place locking rack included with each case, 10 resealable bags of 50 tubes	500	\$141.48
V3802	V4802	Cryogenic Vials, With silicone o-ring, 2mL, sterile, printed graduations and marking area, self-standing star foot, 50-place locking rack included with each case, 10 resealable bags of 50 tubes	500	\$149.04
V3803	V4803	Cryogenic Vials, With silicone o-ring, 3mL, sterile, printed graduations and marking area, self-standing star foot, 50-place locking rack included with each case, 10 resealable bags of 50 tubes	500	\$189.54
V3805	V4805	Cryogenic Vials, With silicone o-ring, 5mL, sterile, printed graduations and marking area, self-standing star foot, 50-place locking rack included with each case, 10 resealable bags of 50 tubes	500	\$213.84

Item	Description	Qty	Price
V5809-A	Cap inserts for cryogenic vials, assorted colors	500	\$18.00
V5809-B	Cap inserts for cryogenic vials, blue	500	\$18.00
V5809-G	Cap inserts for cryogenic vials, green	500	\$18.00
V5809-R	Cap inserts for cryogenic vials, red	500	\$18.00
V5809-W	Cap inserts for cryogenic vials, white	500	\$18.00
V5809-Y	Cap inserts for cryogenic vials, yellow	500	\$18.00
V5920	Locking rack for cryogenic vials, 50-place (5 x 10)	4	\$49.00



# MTC® Erlenmeyer Flasks

- Vented screw caps
- Flat or baffled bottoms
- Gamma sterile (SAL 10<sup>-6</sup>)
- RNase, DNase, pyrogen free
- Individually wrapped



Item	Description	Qty	Price
<b>Flat bottom, disposable, PETG</b>			
F4060-F	125mL Erlenmeyer flask, PETG, flat base	24	\$213.60
F4061-F	250mL Erlenmeyer flask, PETG, flat base	12	\$120.00
F4062-F	500mL Erlenmeyer flask, PETG, flat base	12	\$154.80
F4063-F	1000mL Erlenmeyer flask, PETG, flat base	6	\$117.00
<b>Baffled bottom, autoclavable/reusable, polycarbonate</b>			
F4060-B	125mL Erlenmeyer flask, polycarbonate, baffled base	24	\$239.00
F4061-B	250mL Erlenmeyer flask, polycarbonate, baffled base	12	\$136.00
F4062-B	500mL Erlenmeyer flask, polycarbonate, baffled base	12	\$179.00
F4063-B	1000mL Erlenmeyer flask, polycarbonate, baffled base	6	\$124.40
<b>Baffled bottom, Jumbo flasks, co-polymer</b>			
F4086	2.8L Jumbo flask, styrene/butadiene, baffled base	6	\$270.00
F4087	5.0L Jumbo flask, styrene/butadiene, baffled base	4	\$380.00

# ReadyStrain II™ Sterile Cell Straining Kits

- Includes strainer, 50mL tube & cap
- Opens directly under the hood reduces contamination risks
- Preassembled – ready to use
- Guaranteed sterile – gamma
- Packaged in easy-open “tear-strip” bags



Item	Description	Qty	Price
C5040	ReadyStrain™ 40µm Cell straining kit w/ strainer, 50mL tube and screw-cap, blue, sterile, indiv. wrapped	50	\$166.80
C5070	ReadyStrain™ 70µm Cell straining kit w/ strainer, 50mL tube and screw-cap, white, sterile, indiv. wrapped	50	\$166.80
C5100	ReadyStrain™ 100µm Cell straining kit w/ strainer, 50mL tube and screw-cap, yellow, sterile, indiv. wrapped	50	\$166.80

# SureStrain™ Cell Strainers

- Advanced design – faster flow through
- Compatible with 6-well plates
- Interlocking and stackable
- Sterile; RNase, DNase, pyrogen and cytotoxin free
- Four 1.5mm feet for straining into flat plates and dishes
- Designed for use with 50mL tubes
- Adapter for 15mL tubes included\*



\* Patent Pending

Item	Description	Qty	Price
40CELL	40µm cell strainer, blue, w/ 1 reducing adapter	50	\$109.00
70CELL	70µm cell strainer, white, w/ 1 reducing adapter	50	\$109.00
100CELL	100µm cell strainer, yellow, w/ 1 reducing adapter	50	\$109.00
PEST2	PolyPestle™ double-sided pestle, for use with cell strainers or microcentrifuge tubes, sterile, individually wrapped	50	\$51.00
CELLADAPT	Reducing adapter for use with 15mL, 17x100mm, 12x75mm and 5mL conical tubes, (Patent Pending)	25	\$47.88

# HALO™ Pipettor Single Channel Adjustable Pipettes



## ERGONOMIC DESIGN

Lightweight with stress point reducing contoured edges and soft-ejection, HALO pipettes are designed for comfort.

## ULTRA LIGHT WEIGHT

Up to 40% lighter than the market leaders, HALO™ pipettors are among the lightest ever produced.

## SIMPLE CALIBRATION

In-lab calibration is made possible with the calibration key supplied with every pipette.

## SUPERIOR ACCURACY

With a 4 digit micrometer, and overlapping volume ranges, aspirations of even fractional volumes are achievable with repeatable precision.

## 3 YEAR WARRANTY

Feel confident knowing that your pipette is built to last and covered by a three year manufacturer's warranty.



Receive Custom Laser Engraving with Each HALO Starter Kit

Item	Description	Qty	Price
H6800-2	HALO™ Single Channel Pipettor, lightweight, 0.1-2.0µL	1	\$249.00
H6800-10	HALO™ Single Channel Pipettor, lightweight, 0.5-10µL	1	\$249.00
H6800-20	HALO™ Single Channel Pipettor, lightweight, 2-20µL	1	\$249.00
H6800-50	HALO™ Single Channel Pipettor, lightweight, 5-50µL	1	\$249.00
H6800-100	HALO™ Single Channel Pipettor, lightweight, 10-100µL	1	\$249.00
H6800-200	HALO™ Single Channel Pipettor, lightweight, 20-200µL	1	\$249.00
H6800-1M	HALO™ Single Channel Pipettor, lightweight, 100-1000µL	1	\$249.00
H6800-1MXL	HALO™ Single Channel Pipettor, lightweight, 200-1250µL	1	\$249.00
H6800-5M	HALO™ Single Channel Pipettor, lightweight, 0.5-5mL	1	\$249.00
H6800-10M	HALO™ Single Channel Pipettor, lightweight, 1-10mL	1	\$249.00
H6800-SK2	HALO™ Starter Kit, Includes H6800-10, H6800-20, H6800-200, H6800-1M and Carousel Stand, with FREE Laser Engraving Personalization on Pipettes	1	\$699.00

# PR1MA™ Pipette Controllers with Retractable QuickStand™

- With QuickStand™ retractable stand
- Variable speed control
- Cordless and rechargeable
- 8 hours of continuous operation
- For pipettes from 1mL to 100mL



# PR1MA™ ProPette REACH™ Long-Neck Pipette Controller


- Extended neck for ergonomics
- Perfect for use under the hood
- Allows user to rest elbow on bench
- QuickStand™ retractable kickstand
- 8-hour rechargeable battery life
- Pressure-sensitive speed control



Item	Description	Qty	Price
MPIP-1	PR1MA™ ProPette™ Pipette Controller, w/ QuickStand™	1	\$317.90
MID6950	PR1MA™ ProPette REACH™, Long-Neck Pipette Controller, w/ QuickStand™	1	\$369.50

# Pure • Amp™ PCR Plates


- RNase, DNase, DNA, and PCR inhibitor free
- Compatible with standard equipment
- Ultra-thin walled tube wells
- Alpha-numeric reference grids



Item	Description	Qty	Price
P9601-N	PCR Plates 96 x 0.1ml (Low Profile/Fast) non-skirted, 50/pk	50	\$173.68
P9601-NW	PCR Plates 96 x 0.1ml (Low Profile/Fast) non-skirted, WHITE, 50/pk	50	\$196.56
P9601-S	PCR Plates 96 x 0.1ml (Low Profile/Fast) Semi Skirted, 50/pk	50	\$190.91
P9601-SW	PCR Plates, 96 x 0.1ml (Low Profile/Fast) Semi Skirted, WHITE, 50/pk	50	\$216.54
P9601-SRR	PCR Plates 96 x 0.1ml Raised Rim (Low Profile/Fast) Semi Skirted, 50/pk	50	\$178.10
P9602-N	PCR Plates, Standard 96 well x 0.2ml, non-skirted, 50/pk	50	\$128.13
P9602-NW	PCR Plates, Standard 96 well x 0.2ml, non-skirted, WHITE, 50/pk	50	\$161.91
P9602-S	PCR Plates, Semi Skirted, 96 x 0.2ml, 50/pk	50	\$160.16
P9602-SRR	PCR Plates, Semi Skirted with Raised Rim (ABI), 96 x 0.2ml, 50/pk	50	\$160.16

# Pure • Amp™ Pre-cut Sealing Membranes

- RNase, DNase, DNA, and PCR inhibitor free
- Evaporation rates ≤ 5%
- Inner packs of 10 membranes



Item	Description	Qty	Price
P1001-A	Sealing film, Aluminum membrane	100	\$82.68
P1001-PCR	Sealing film, PCR 96 Well Plate Sealing Membrane, 10 bags of 10 seals	100	\$82.68
P1001-Q	Sealing film, qPCR Optical (sticky adhesive film) Bio-Rad Type, 10 bags of 10 seals	100	\$121.18
P1001-QPS	Sealing film, qPCR Optical (pressure-sensitive adhesive) ABI Type, 10 bags of 10 seals	100	\$165.80
P1001-E	Sealing film, ELISA plates, 10 bags of 10 seals	100	\$82.68

# Standard Solution Reservoirs

- Graduated inner walls
- Micro-trough™ U-shaped bottoms
- Pouring spouts on all 4 corners
- RNase, DNase, and pyrogen free
- One pack of 5 Protective Lids (P8025-LID) included with standard 25ml reservoirs



Item	Description	Qty	Price
P8025	Reservoir, 25ml non-sterile, bulk pack, includes 5 lids	100	\$58.24
P8025-1S	Reservoir, 25ml sterile, Individually wrapped, includes 5 lids	100	\$122.08
P8025-5S	Reservoir, 25ml sterile, 5/ sterile bag, includes 5 lids	200	\$122.08
P8025-LID	Lids for 25ml solution reservoirs, sterile, 5 per bag	200	\$84.53
P8050	Reservoir, 50ml non-sterile, bulk pack	100	\$58.24
P8050-1S	Reservoir, 50ml sterile, Individually wrapped	100	\$122.08
P8050-5S	Reservoir, 50ml sterile, 5/ sterile bag	200	\$122.08
P8100	Reservoir, 100ml non-sterile, bulk pack	100	\$58.24
P8100-1S	Reservoir, 100ml sterile, Individually wrapped	100	\$122.08
P8100-5S	Reservoir, 100ml sterile, 5/sterile bag	200	\$122.08
P8051	Reservoir, divided, 50ml non-sterile, bulk pack	100	\$58.24
P8051-1S	Reservoir, divided, 50ml sterile, Individually wrapped	100	\$122.08
P8051-5S	Reservoir, divided, 50ml sterile, 5/sterile bag	200	\$122.08

# Pure • Amp™ PCR Tubes, Strips & Caps

- RNase, DNase, DNA, and PCR inhibitor free
- 0.2mL or Low Profile 0.1 mL options
- Ultra-thin, uniform wall thickness
- White and frosted tubes for enhanced qPCR results



Item	Description	Qty	Price
P3010-A	PCR 8-Strips. 0.2mL with Attached Individual Flat Caps	120	\$82.68
P3010-S	PCR 8-Strips, 0.2mL with Separate Flat Cap Strips	120	\$82.68
P3020	PCR 8-Strips, 0.2mL with Separate Dome Cap Strips	120	\$82.68
P3030	0.2mL PCR tubes, individual, with flat caps	1000	\$47.84
P3040	0.2mL PCR tubes, individual, with domed caps	1000	\$47.84
P3801-Q	0.1 mL qPCR 8-Strip (With SEPARATE Optical Strip Caps), Natural/Clear	120	\$142.32
P3801-QF	0.1 mL qPCR 8-Strip (With SEPARATE Optical Strip Caps) FROSTED	120	\$150.54
P3801-QW	0.1 mL qPCR 8-Strip (With SEPARATE Optical Strip Caps), WHITE	120	\$150.54
P3802-Q	0.2 mL qPCR 8-Strip (With SEPARATE Optical Strip Caps), Natural/Clear	120	\$142.32
P3802-QW	0.2 mL qPCR 8-Strip (With SEPARATE Optical Strip Caps), WHITE	120	\$150.54
PR-PCRTOOL	Tool for opening/closing PCR strips	2	\$86.80

# ASPIR-8™ Solution Reservoirs

- Specially designed for 8-Channel pipettes
- 75mm inner chamber
- Conserve valuable reagents



Item	Description	Qty	Price
P7005	ASPIR-8™ Reservoir, 5mL, non-sterile, bulk pack	300	\$111.28
P7005-1S	ASPIR-8™ Reservoir, 5mL, sterile, Individually wrapped	100	\$111.28
P7005-5S	ASPIR-8™ Reservoir, 5mL, sterile, 5/ sterile bag	200	\$111.28
P8010	ASPIR-8™ Reservoir, 10mL, non-sterile, bulk pack	300	\$111.28
P8010-1S	ASPIR-8™ Reservoir, 10mL, sterile, Individually wrapped	100	\$111.28
P8010-5S	ASPIR-8™ Reservoir, 10mL, sterile, 5/ sterile bag	200	\$111.28
H66725	ASPIR-8™ Reservoir, 25mL, non-sterile, bulk pack	100	\$55.64
H66725-1S	ASPIR-8™ Reservoir, 25mL, sterile, Individually wrapped	100	\$111.28
H66725-5S	ASPIR-8™ Reservoir, 25mL sterile, 5/ sterile bag	200	\$111.28

# BenchBin™

## Benchtop Biohazard Bin

- Safely dispose of tips, tubes, gloves, etc.
- Secures and covers hazardous waste
- Easy push open/close spring lid
- Compact footprint – 7 x 5 x 7in.



Item	Description	Dimensions	Qty	Price
A8401B	BenchBin™ Benchtop Biohazard Bin, 1.5L	7 x 5 x 7in. (178 x 127 x 178mm)	1	\$39.80
A9004R	Red Biohazard Bags w/ black printed markings	8.5 x 11in. (216 x 279mm)	100	\$31.00
A8401-SK	BenchBin™ Starter Kit, Includes 1 bin and 400 autoclave bags (8.5 x 11in.)		1	\$115.00

# The BowTie™

## Collection

### Biohazard Bags & Bin

- Drawstring for instant closure
- 12-16 gallon capacity, 25 x 35in (45-60 liter capacity, 64 x 89cm)
- Color changing sterilization indicator
- Bags fit perfectly in Biohazard Bin
- Bin has a hands-free foot pedal operated lid



Item	Description	Qty	Price
A8001R	BowTie™ Biohazard bags, with drawstring, HDPE, 12-16 gallon capacity (45-60 liters), 25 x 35in. (64 x 89cm), large marking area and sterilization indicator	100	\$159.43
A8000B	BowTie™ Biohazard Bin, with hands-free foot pedal, PP, 12-14 gallon capacity (45-52 liters), 10 x 14 x 19in.. (25 x 35 x 48cm, W x L x H), attached lid	1	\$74.79

# SureSeal S™ Sterile MicroTubes

- Sterile comes standard – by electron beam radiation
- Self-Standing bag – with gusseted bottom
- Stop-Pops™ supplied – 10 cap locks in every pack
- All SureSeal S tubes are RNase, DNase and pyrogen free.



Stop-Pops™



Self-Standing Bag

Item	Description	Qty	Price
MID15C	Microtube w/ cap, 1.5mL, clear, sterile, w/ self-standing bag & Stop-Pops™	500	\$17.01
MID15C-OB	Microtube w/ cap, 1.5mL, opaque black, sterile, w/ self-standing bag & Stop-Pops™	500	\$17.40
MID15C-50	Microtube w/ cap, 1.5mL, clear, sterile, w/ self-standing bag, 10 packs of 50	500	\$20.26
MID15C-AST	Microtube w/ cap, 1.5mL, assorted colors, (B, R, G, O, P, Y), sterile, w/ self-standing bag & Stop-Pops™	500	\$19.37
MID15C-50AST	Microtube w/ cap, 1.5mL, assorted colors, sterile, w/ self-standing bag, 10 single color bags of 50	500	\$20.26
MID2C	Microtube w/ cap, 2.0mL, clear, sterile, w/ self-standing bag & Stop-Pops™	500	\$19.08
MID05C	Microtube w/ cap, 0.5mL, clear, sterile, w/ self-standing bag	500	\$19.08
MID05C-AST	Microtube w/ cap, 0.5mL, assorted (B, R, G, P, Y), sterile, w/ self-standing bag	500	\$20.28

# ClearSeal™

## Screw Cap Microtubes

- Ultra-clear polypropylene
- Packaged in resealable bags of 50
- Leak-proof O-ring seal
- Sterile, RNase, DNase & Pyrogen free
- White, printed graduations



0.5 mL Conical



0.5 mL Self-Standing



1.5 mL Conical



1.5 mL Self-Standing



2.0 mL Conical



2.0 mL Self-Standing



2.0 mL Self-Standing Loop Cap

Item	Description	Qty	Price
C3205-C	ClearSeal™ screw-cap microtubes, 0.5mL, w/ O-ring & caps, sterile, non-graduated, conical bottom, 20 resealable bags of 50 tubes	1000	\$245.57
C3205-S	ClearSeal™ screw-cap microtubes, 0.5mL, w/ O-ring & caps, sterile, non-graduated, self-standing, 20 resealable bags of 50 tubes	1000	\$245.57
C3215-CG	ClearSeal™ screw-cap microtubes, 1.5mL w/ O-ring & caps, sterile, non-graduated, conical bottom, 20 resealable bags of 50 tubes	1000	\$245.57
C3215-SG	ClearSeal™ screw-cap microtubes, 1.5mL, w/ O-ring & caps, sterile, printed graduations, self-standing, 20 resealable bags of 50 tubes	1000	\$245.57
C3220-CG	ClearSeal™ screw-cap microtubes, 2.0mL, w/ O-ring & caps, sterile, printed graduations, conical bottom, 20 resealable bags of 50 tubes	1000	\$245.57
C3220-SG	ClearSeal™ screw-cap microtubes, 2.0mL, w/ O-ring & caps, sterile, printed graduations, self-standing, 20 resealable bags of 50 tubes	1000	\$245.57
C3170	ClearSeal™ screw-cap microtubes, 2.0mL, w/ O-ring & caps, sterile, generic, non-graduated, colorless caps, self-standing, 20 resealable bags of 50 tubes	1000	\$201.16
C3220-SGL	ClearSeal™ screw-cap microtubes, 2.0mL, w/ O-ring & loop-cap, sterile, printed graduations, self-standing, 20 resealable bags of 50 tubes	1000	\$256.27
R1081-A	BetterBox™ Storage Box with hinged lid, 81 x 1.5mL to 2.0mL Microtubes pack of 5, assorted colors	5	\$42.27

This promo cannot be combined with other promotions or discounts. All terms apply. Promo valid until 06/30/2023. Visit [www.midsci.com](http://www.midsci.com) for more information.



ORDER TODAY! [midsci.com](http://midsci.com) • [custserv@midsci.com](mailto:custserv@midsci.com) • 800.227.9997 • f: 636.225.9998