TECHNICAL INFORMATION



MICEBIO GRYOGERIG VIALS



- Temps as low as -196°C
- Self-standing star foot
- Leakproof silicone O-ring
- Resealable bags of 50
- Gamma sterilized
- Free locking rack
- RNase, DNase, pyrogen, & endotoxin free
- Autoclavable



SPECIFICATIONS

Material LN₂ grade (tube/cap): polypropylene

Material Silicone (o-ring/washer):

Temp Range: -196°C to +121°C*

Vial Diameter: 12.5mm

*LN2 gas phase only; not for liquid phase



CRYOGENIC VIALS									
ITEM #	DESCRIPTION	CAP THREAD	HEIGHT (WITH CAP)	VOL.	QTY	LIST PRICE			
V3801	Cryo Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	External	40 mm	1 mL	500/cs	\$131.00			
V3802	Cryo Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	External	47.3 mm	2 mL	500/cs	\$138.00			
V3803	Cryo Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	External	66.6 mm	3 mL	500/cs	\$175.50			
V3805	Cryo Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	External	90.7 mm	5 mL	500/cs	\$198.00			
V4801	Cryo Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	Internal	40.2 mm	1 mL	500/cs	\$131.00			
V4802	Cryo Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	Internal	47.5 mm	2 mL	500/cs	\$138.00			
V4803	Cryo Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	Internal	60.5 mm	3 mL	500/cs	\$175.50			
V4805	Cryo Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	Internal	91 mm	5 mL	500/cs	\$198.00			
V5809-A*	* Color Cap inserts for Cryogenic Vials, assorted colors (white, blue, green, red, yellow), 500/pk				500/pk	\$18.00			
V5920	Locking Rack for 50 "star foot" vials, 210 X 100 X 25mm				4/pk	\$49.00			

ACCESSORIES								
ITEM #	DESCRIPTION	QTY	LIST PRICE					
R1060	Freezer Box, slip on lid, polypropylene, 81 X 2 ml, 130 X 30 X 45 mm	5/pk	\$42.00					
R2281	Freezer Box, hinged lid, polypropylene, 81 X 2 ml, 130 X 30 X 53 mm	5/pk	\$59.88					
R2200	Freezer Box, hinged lid, polypropylene, 100 X 2 ml, 130 X 30 X 53 mm	5/pk	\$59.88					
R2781-A*	FlipTop™ Freezer Box, hinged lid, cardboard, 81 X 2mL, 133 X 133 X 53mm	5/pk	\$17.70					
R2700-A*	FlipTop™ Freezer Box, hinged lid, cardboard, 100 X 2mL, 133 X 133 X 53mm	5/pk	\$18.42					
R3781	FlipTop™ Freezer Box, hinged lid, cardboard, 81 X 3-4mL, 133 X 133 X 75mm	50/cs	\$148.73					
R3700	FlipTop™ Freezer Box, hinged lid, cardboard, 100 X 3-4mL, 133 X 133 X 75mm	50/cs	\$148.73					



