

# Macro Volume Tip Compatibility Chart

Tips shown in actual size.

## 5 mL Gilson / Universal Fit

Item #	Description	Color	Packaging	Eppendorf	Socorex	Gilson	Rainin Classic	Sartorius (Biohit)	Finnpipette	Brand / Branchtech
4410-00	Bulk	Clear	250/unit, 10 units/case							
4410-0F	Bulk, Filtered	Clear	100/unit, 10 units/case							
4411-00	Racked	Clear	50 tips/unit, 10 units/case			•	•			
4411-50	Racked, Sterile	Clear	50 tips/unit, 10 units/case							
4411-SF	Racked, Sterile, Filtered	Clear	50 tips/unit, 10 units/case							

## 5 mL Eppendorf Fit

Item #	Description	Color	Packaging	Eppendorf	Socorex	Gilson	Rainin Classic	Sartorius (Biohit)	Finnpipette	Brand / Branchtech
4420-00	Bulk	Clear	250/unit, 10 units/case							
4420-0F	Bulk, Filtered	Clear	100/unit, 10 units/case							
4421-00	Racked	Clear	50 tips/unit, 10 units/case	•	•			•*		
4421-50	Racked, Sterile	Clear	50 tips/unit, 10 units/case							
4421-SF	Racked, Sterile, Filtered	Clear	50 tips/unit, 10 units/case							

## 5 mL Thermo Finn / Brand Fit

Item #	Description	Color	Packaging	Eppendorf	Socorex	Gilson	Rainin Classic	Sartorius (Biohit)	Finnpipette	Brand / Branchtech
4440-00	Bulk	Clear	5 bags of 100/unit, 2 units/cs							
4441-00	Racked	Clear	50 tips/unit, 10 units/case						•	•
4441-50	Racked, Sterile	Clear	50 tips/unit, 10 units/case							
4441-SF	Racked, Sterile, Filtered	Clear	50 tips/unit, 10 units/case							

## 10 mL Thermo Finn / Socorex Fit

Item #	Description	Color	Packaging	Eppendorf	Socorex	Gilson	Rainin Classic	Sartorius (Biohit)	Finnpipette	Brand / Branchtech
4520-01	Bulk	Blue	200/unit, 5 units/case							
4521-01	Racked	Blue	25 tips/rack, 10 racks/case		•				•	
4521-S1	Racked, Sterile	Blue	25 tips/rack, 10 racks/case							
4521-SF	Racked, Sterile, Filtered	Blue	25 tips/rack, 10 racks/case							

## 10 mL Gilson / Eppendorf Fit

Item #	Description	Color	Packaging	Eppendorf	Socorex	Gilson	Rainin Classic	Sartorius (Biohit)	Finnpipette	Brand / Branchtech
4530-00	Bulk	Clear	200/unit, 5 units/case							
4530-0F	Bulk, Filtered	Clear	50/unit, 10 units/case							
4531-00	Racked	Clear	25 tips/rack, 10 racks/case	•		•	•	•*		•
4531-50	Racked, Sterile	Clear	25 tips/rack, 10 racks/case							
4531-SF	Racked, Sterile, Filtered	Clear	25 tips/rack, 10 racks/case							

\*Does not fit Proline pipettors, but does fit all other Biohit pipettors, including Proline Plus.

