

## Non-sterile Powder Free Nitrile Examination Gloves

	YTY INDUSTRY (MA	NJUNG) SDN BI	HD	
Spec No. YTY1004/03 – 5ASTM	PRODUCT SPECIFICATION		Eff	Revision: 3 fective: 22/02/2017
Product Information:	•	<u>.</u>		
Trade Name	Non-sterile Powder Free Nitrile Examination Gloves			
Intended Use	A patient examination glove is a disposable device intended for medical purpose that is worn on the examiner's hand or finger to prevent contamination between patient and			
Product Conformance	examiner ASTM D6319 Current Revision and Customer Requirement			
Material	Nitrile (Carboxylated Butadiene Acrylonitrile)			
Cuff Finishing	Beaded			
Color	Blue			
Туре	Ambidextrous			
Surface finish	Finger textured			
Product Specification:	1 0			
Dimension	Size	Palm Width (mm)		Length (mm)
	X-Small	75 ± 5mm		240min
	Small	85 ± 5mm		240min
	Medium	95 ± 5mm		240min
	Large	106 ± 5mm		240min
	Extra Large	115 ± 5mm		240min
	XX Large	125 ± 5mm		240min
Thickness	701 Lange	Single wall thickness (mm)		
		Specification		
	Cuff		0.05-0.08	
	Palm			-0.09
	Finger			-0.14
				imin
	Bodding	Un-aged Aged		
Physical Properties (6mm die cut)		Specification		Specification
,	Tensile Strength (MPa)	14min		14min
	Ultimate Elongation (%)	500min		400min
Residual Powder		2.0 mg/glo	ve maximum	
Quality Specification  * AND – Accept No Defects	Inspection Criteria			Sampling Plan/AQL
	Dimension			S2 – AQL 4.0
	Physical Properties			S2 – AQL 4.0
	Water Tight Test			G1 – AQL 1.0
				(Single Normal Sampling)
	Visual Inspection – Major			G1 – AQL 2.5
	Visual Inspection – Minor			G1 – AQL 4.0
	Residual Powder			n = 5 (average)
	Packaging Labeling			n = 5 (AND)
Destruction O. C.	Glove Physical Counting			n = 5 (AND)
Packaging Configuration and Labeling				
Packing	Gloves per Inner Dispenser (by weight) is per customer's requirement (e.g 100/200/250 10 Dispensers per Outer Carton			
Labeling	The dispensers and outer cartons shall have the following information/prints: Product description and Brand including size and net quantity of gloves. Country of Origin of Manufacturer and name and address of distributor. Single Use Only and/or the Symbol for single use and Non-Sterile.			
Lot Numbering System	Per the customer's requirement or in house as per packing SOP.  Manufacturing and/or expiry date where applicable.			
Authorized by:	<u>.                                      </u>			Code:
Chow Choong Keong QA Manager				N35TNB [8612F (C)]
WA Ivialiayei				