


## PR1MA-Shell Rigid PCR Plates

- Rigid polycarbonate frame to reduce warping commonly associated with polypropylene PCR plates.
- Ultra-thin polypropylene wells provide optimal thermal transfer and higher reaction efficiencies.
- White well options enhance the fluorescence signal during qPCR.
- Broad thermal cycler compatibility to fit most commercially-available instruments.
- Raised rim around each well for improved sealing.
- High contrast printed alphanumeric.
- ANSI/SBS footprint and well spacing make them perfectly suitable for liquid handling instruments.
- Certified DNase, RNase, and pyrogen-free.

Online PCR  
Consumables  
Selector



Item #	# of Wells	Volume	Description	Frame Color	Well Color	Packaging	List Price
<b>PR-PCR3742-C</b>	96	0.1 mL	Full-Skirt	Clear*	Clear	10 plates/pk	\$56.78
<b>PR-PCR3742-CW</b>	96	0.1 mL	Full-Skirt	Clear	White	10 plates/pk	\$56.78
<b>PR-PCR3750-C</b>	96	0.2 mL	Semi-Skirt	Clear*	Clear	10 plates/pk	\$56.78
<b>PR-PCR3735-C</b>	384	30 µL	Full-Skirt, Single Notch	Clear*	Clear	10 plates/pk	\$65.30
<b>PR-PCR3735-CW</b>	384	30 µL	Full-Skirt, Single Notch	Clear	White	10 plates/pk	\$83.77
<b>PR-PCR3730-C</b>	384	30 µL	Dual-Notch	Clear	Clear	10 plates/pk	\$79.98

\* Available in individual colors (blue, green, red, violet, yellow) or in an assortment of the five colors.

## PR1MA Polypropylene PCR Plates

- Broad thermal cycler compatibility to fit most commercially-available instruments.
- Raised rim around each well for improved sealing.
- Optically clear with printed black alphanumeric.
- Ultra-thin polypropylene wells provide optimal thermal transfer and higher reaction efficiencies.
- Certified DNase, RNase, and pyrogen-free.




Item #	# of Wells	Volume	Description	Color	Notch	Packaging	List Price
<b>PR-PCR1296</b>	96	0.2 mL	No-Skirt	Clear	A12	10 plates/pk	\$24.80
<b>PR-PCR1196</b>	96	0.1 mL	No-Skirt, Low-Profile	Clear	H12	10 plates/pk	\$49.61
<b>PR-PCR2296S</b>	96	0.2 mL	Semi-Skirt	Clear	A12	10 plates/pk	\$24.80
<b>PR-PCR2296</b>	96	0.2 mL	Semi-Skirt, ABI®-Type	Clear	A12	10 plates/pk	\$24.80
<b>PR-PCR2296R</b>	96	0.2 mL	Semi-Skirt, ABI®-Type, Raised Rim	Clear	A12	10 plates/pk	\$24.80
<b>PR-PCR2196F</b>	96	0.1 mL	Semi-Skirt, Low Profile, Raised Sides, ABI Fast®-Type	Clear	A1	10 plates/pk	\$24.80
<b>PR-PCR2196LC-W</b>	96	0.1 mL	Semi-Skirt, Low Profile, LightCycler® Compatible	White	H12	10 plates/pk	\$29.85
<b>PR-PCR3196</b>	96	0.1 mL	Full-Skirt, Low-Profile	Clear	A12	10 plates/pk	\$24.80
<b>PR-PCR3196-W</b>	96	0.1 mL	Full-Skirt, Low-Profile, Bio-Rad® Compatible	White	A12	10 plates/pk	\$29.85
<b>PR-PCR33841</b>	384	30 µL	Full-Skirt, ABI®-Type	Clear	A24	10 plates/pk	\$36.20
<b>PR-PCR33842</b>	384	30 µL	Full-Skirt, Roche®-Type	Clear	A24, P24	10 plates/pk	\$35.67



## PCR Strip Tubes



**Combo Pack**  
PR-PCR28CF



Rigid Strips feature robust, **3-point interwell linkages** making strips more sturdy and decreasing the chance of spillage.



**Break-Away Strips** feature single interwell linkages, making strips easy to split into smaller sections.



**Optically-clear caps** ensure that the fluorescence signal is passed through uninhibited.

Item #	Type	Volume	Description	Packaging	List Price
PR-PCR28ACF	Rigid	0.2 mL	8-Strip Tubes with Individually Attached Hinge Caps	120 strips/pk	\$68.86
PR-PCR18ACF	Rigid	0.1 mL	8-Strip Tubes with Individually Attached Hinge Caps, Low-Profile	120 strips/pk	\$70.84
PR-PCR18ACF-W	Rigid	0.1 mL	8-Strip Tubes with Individually Attached Hinge Caps, Low-Profile, White	120 strips/pk	\$74.14
PR-PCR28ACD	Break-Away	0.2 mL	8-Strip Tubes with Individually Attached Dome Caps	120 strips/pk	\$68.86
PR-PCR28ACFF	Break-Away	0.2 mL	8-Strip Tubes with Individually Attached Flat Caps	120 strips/pk	\$68.86
PR-PCR28CD	Combo Pack	0.2 mL	8-Strip Tubes and Dome Caps	125 strips/pk	\$76.82
PR-PCR28CF	Combo Pack	0.2 mL	8-Strip Tubes and Flat Caps	125 strips/pk	\$76.82
PR-PCR18CF	Combo Pack	0.1 mL	8-Strip Tubes and Flat Caps, Low-Profile	125 strips/pk	\$76.82

## Individual PCR Tubes

Item #	Volume	Description	Packaging	List Price
PR-PCR2D	0.2 mL	Dome Cap PCR Tube	1,000/pk	\$45.33
PR-PCR2F	0.2 mL	Optically-Clear, Flat Cap PCR Tube	1,000/pk	\$38.80
PR-PCR2FF	0.2 mL	Frosted, Flat Cap PCR Tube	1,000/pk	\$38.80



**Online PCR Consumables Selector**

### Your MIDSCI Sales Representative:

Name: \_\_\_\_\_

Email: \_\_\_\_\_

Phone: \_\_\_\_\_