

# SPEC OF AT500M-9

STP-JX-Tip-AT500M-9-A1

GEB Robotic Tips for Tecan MCA®. Provides reservoirs with the High Throughput Screening (HTS) samples & reagents storage. Used in proteomics, drug development, and genomics.

## Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (4.7*3.5)mm	Packaging Spec.
AT500M-9-N	500ul	Natural				96 tips/rack,
AT500M-9-N-LB	500ul	Natural	Low Retention			
AT500M-9-NS	500ul	Natural		Sterilized		5 packs / case.
AT500M-9-NS-LB	500ul	Natural	Low Retention	Sterilized		50 racks/case, 4800tips/case.
AF500M-9-NS	500ul	Natural		Sterilized	Filtered	
AF500M-9-NS-LB	500ul	Natural	Low Retention	Sterilized	Filtered	
Case dimension L*W*H: 500*285*405mm			Gross Weight : 10.5 kg			

- ◇ 500µL Natural Tips
- ◇ 96 Tip/Hanging Tip Rack format
- ◇ Polypropylene
- ◇ Certified DNase/RNase free
- ◇ Nonpyrogenic
- ◇ Presterile/Nonsterile
- ◇ Filtered (Filtered With aerosol-resistant filters)/NonFiltered
- ◇ Endotoxin-Safe meets USP/FDA standards
- ◇ Comply With USP 3127 (Heavy Metal Contamination) standards
- ◇ Compatibility: Tecan MCA® workstations.

Drawing: (Diameter tolerance:±0.03mm , Height tolerance:±0.25mm)

