



HERMLE
Benchmark



Part# Z496-K

Part# Z496-K-UC

Hermle-Benchmark Z496-K Series

**Compatible with 450ml and 600ml bioreactor bottles
3750 rpm / 2400 xg**



Part# Z496-1000-BR
400-600ml bioreactor rotor and buckets complete set

HERMLE Benchmark



Part# Z446



Part# Z446-K

Hermle-Benchmark Z446 Series

Compatible with 450ml bioreactor bottles
4000 rpm / 3434 xg



Part# Z446-750 + Z446-750-A450 (x2)
4 x 750ml rotor and buckets with 450ml bioreactor buckets