

## SealPlate<sup>®</sup> ColorTab™

For sample identification and plate labeling

## **Application:**

General laboratory use, short-term storage, incubation and ELISA


## **Temperature Range:**

-40°C to +120°C

CATALOG #	DESCRIPTION	QUANTITY
SP-IDA-10	SealPlate ColorTab, Assorted	50
SP-ID"X"-100*	SealPlate ColorTab, Color	100

Each end-tab features a colored dot and white writing area. An inscribed color-coded end-tab can be removed and pasted into a laboratory notebook while the other remains on the plate for identification. The film sheet can be removed from the plate without disturbing the identification label due to the finely perforated end-tabs. Choose from 5 different individual colors or the assortment pack.

- Provides a secure seal and is easy to remove
- Rigid polyester film features good clarity and easy handling
- Features two end-tabs with perforations



"For 3 times faster sealing, check out our SealMate films on a roll with dispenser"

