

# COMFORT AND DEXTERITY FOR EVERYDAY HAND PROTECTION.







## TF-95LG

POWDER-FREE NITRILE EXAMINATION GLOVES  
WITH ENERGEL REJUVINATING TECHNOLOGY

**FEATURING AN ENERGEL INNER COATING,  
THESE GLOVES SOOTHE AND MOISTURIZE DRY  
HANDS WITHOUT SACRIFICING DEXTERITY**

- EnerGel® Hand Conditioning Technology clinically proven to hydrate and soothe skin
- Exceeds 1.5 AQL laboratory standards
- Fully textured surface for excellent wet and dry grip
- Touch-screen compatible
- Tear-resistant beaded cuff eases donning
- SHA Dermatologically Approved


TF					TF-95LG	
WEIGHT	LENGTH	THICKNESS	FEATURES	COLOR		
 5.0g (large size)	 9.5"	 3.6 mil	EnerGel®	 LIGHT GREEN		



Optimal  
thickness  
for durability  
and tactile  
sensitivity

Fully textured  
surface for wet  
and dry grip

Tear-resistant  
beaded cuff for  
easy donning



**TF-95LG**

**POWDER-FREE NITRILE EXAMINATION GLOVES  
WITH ENERGEL REJUVINATING TECHNOLOGY**

**PRODUCT DESCRIPTION****SPECIFICATIONS**

<b>PRODUCT</b>	TrueForm®
<b>GRADE</b>	Powder-Free Exam
<b>MATERIAL</b>	Nitrile
<b>THICKNESS</b>	<b>FINGER</b> 4.0mil 0.10mm
	<b>PALM</b> 3.6mil 0.09mm
	<b>CUFF</b> 3.0mil 0.08mm
<b>COLOR</b>	Light Green
<b>SURFACE</b>	<b>OUTER</b> Fully Textured
	<b>INNER</b> EnerGel®
<b>LENGTH</b>	9.5" 242mm
<b>WEIGHT</b>	5.0g
<b>FIT TYPE</b>	Ambidextrous
<b>CUFF TYPE</b>	Beaded
<b>WATERTIGHT AQL</b>	1.5
<b>ORIGIN COUNTRY</b>	China

**PHYSICAL PROPERTIES**

	PRE-AGING	POST-AGING
<b>TENSILE STRENGTH:</b>	22MPA	20 MPA
<b>ELONGATION:</b>	560%	550%

**INTEGRATED TECHNOLOGIES**

**ENERGEL** Hand Conditioning Technology hydrates and soothes the skin with purified aloe and has been clinically proven to prevent many skin ailments, like dermatitis. It has a healing effect that improves skin health.

Clean Washed

100% Air Inspected



\* Specifications based on size Large glove

**STANDARDS & CERTIFICATIONS**

EN ISO 21420:2020  
EN 455-2000-2015  
ASTM D6319-19  
ASTM F1671-13  
ASTM F739-12  
Reach EU 1907/2006  
EU 2016/425

CAT III Complex Design  
UL NFPA 1999-2018  
FDA 21 CFR 177.2600



Please visit [swsafety.com/symbols](http://swsafety.com/symbols) for information on standards used on this document.

**ORDERING**

ITEM NO.	SIZE	PALM WIDTH	PACKAGING
N128401	XS 5½ – 6	3.0" 75mm	100 gloves per box 10 boxes per case
N128402	S 6½ – 7	3.4" 86mm	100 gloves per box 10 boxes per case
N128403	M 7½ – 8	3.7" 95mm	100 gloves per box 10 boxes per case
N128404	L 8½ – 9	4.2" 107mm	100 gloves per box 10 boxes per case
N128405	XL 9½ – 10	4.6" 117mm	100 gloves per box 10 boxes per case

**APPLICATIONS****INDUSTRIAL**

Assembly  
Liquid Handling

**MEDICAL**

Chemical Handling  
Clinical  
First Responder EMS  
Handling and Packaging Class 1  
Medical Devices  
Laboratory  
Light Chemical Mixing  
Pharma Production

Keep out of sunlight. Store in a cool, dry place.  
Keep away from sources of ozone or ignition.



R\_N12840X1121

