



ThermalSeal® 2x8 MiniStrips™

For sealing one or two rows of a 96-well plate

Application:

Short-term storage


Temperature Range:

-40°C to +120°C

CATALOG #	DESCRIPTION	QUANTITY
TS-2X8-50	ThermalSeal MiniStrips	200

ThermalSeal 2x8 MiniStrips are less expensive than full-size sheets when using only a partial plate. Designed for use with standard 96- well plates where rows must be selectively protected or accessed.

- Certified DNase- & RNase- free
- Each sheet contains 4 strips on a continuous liner
- Apply individually or two-at-a-time
- Can be used as a secondary film to seal previously pierced wells or to prevent accidental pipetting



Single 2x8 Strip

"For a 1x8 PCR FilmStrip featuring an extra strong adhesive, check out EZcap 1x8 FilmStrips"

