



MIDSCI PCR Plate-Thermal Cycler Fit Chart

Thermal Cycler	Non-Skirted																Low Profile	
	AVRT-N	B70501	AVT2619	B50501	AB70651	AVRT1	B70651	AVT7437	B50651	AVT3887	B17480	AB77501	AVT3890	AVS906	AVT2844	B70515L		B71559L
ABI	2400																	2400
	2700																	2700
	2720																	2720
	5700, 7000																	5700, 7000
	7300, 7500, 7700																	7300, 7500, 7700
	7500 Fast																	7500 Fast
	7900																	7900
	7900HT																	7900HT
	9600, 9700																	9600, 9700
	9800 Fast																	9800 Fast
	StepOne Plus																	StepOne Plus
	Veriti 96																	Veriti 96
	Via7																	Via7
BioRad	GeneCycler																	GeneCycler
	Cycler																	Cycler
	Cycler IQ																	Cycler IQ
	iQ5																	iQ5
	MyCycler																	MyCycler
	MyiQ																	MyiQ
	Ptc-100, 200, 225																	Ptc-100, 200, 225
	Chrom4, CFX96																	Chrom4, CFX96
	Opticon, Opticon2																	Opticon, Opticon2
	C-1000, S-1000																	C-1000, S-1000
	T100																	T100
Benchmark Scientific	BSTC-96																	BSTC-96
	BSTC-384																	BSTC-384
DyNA	DCS, DCR																	DCS, DCR
Eppendorf	MasterCycler																	MasterCycler
	MasterCycler Gradient																	MasterCycler Gradient
	MasterCycler ep Gradient/Gradient S, Pro																	MasterCycler ep Gradient/Gradient S, Pro
	MasterCycler ep Realplex																	MasterCycler ep Realplex
Roche	Light Cycler 480																	Light Cycler 480
Labnet International	MultiGene Optimax																	MultiGene Optimax
Stratagene	Robo-Cycler																	Robo-Cycler
Bibby Scientific	MX3000P, 3005P, 4000																	MX3000P, 3005P, 4000
	3Prime, Prime, Prime-Elite																	3Prime, Prime, Prime-Elite
	PrimeQ																	PrimeQ

MIDSCI PCR Plate- Sealing Film/Strip Cap Recommended Pairings

STERILE	PERCABLE	PCR FILM	Non-Skirted																Low Profile
			AVRT-N	B70501	AVT2619	B50501	AB70651	AVRT1	B70651	AVT7437	B50651	AVT3887	B17480	AB77501	AVT3890	AVS906	AVT2844	B70515L	
		T52-100																	T52-100
		100-THER-PLT																	100-THER-PLT
		STR-THER-PLT																	STR-THER-PLT
		AF-100																	AF-100
		AFS-25																	AFS-25
		E-384-100																	E-384-100
		TSA-100																	TSA-100
		F-96-100																	F-96-100
		TS-2X8-50																	TS-2X8-50
		PCR FILM																	PCR FILM
		T55-RTQ-100																	T55-RTQ-100
		TS-RT2-100																	TS-RT2-100
		TS-RT2RR-100																	TS-RT2RR-100

12 STRIP CAP	8 STRIP CAP	FLAT CAP	DOME CAP	STRIP CAPS	Non-Skirted																Low Profile				
					AVRT-N	B70501	AVT2619	B50501	AB70651	AVRT1	B70651	AVT7437	B50651	AVT3887	B17480	AB77501	AVT3890	AVS906	AVT2844	B70515L	B71559L	PCR Natural Plates	PCR Frosted Plates	PCR White Plates	
					B79701																				B79701
					C79701																				C79701
					B57801																				B57801
					AVSTC-N																				AVSTC-N
					AVSTC-12																				AVSTC-12
					AVT3799																				AVT3799
					AVT3876																				AVT3799

MIDSCI Strip Tube-Strip Cap Fit Chart

12 Strip Cap	8 Strip Cap	Flat Cap	Dome Cap	Strip Tubes										
				AVRT-N	B70501	AVT2619	B50501	AB70651	AVRT1	B70651	AVT7437			
				B79701										
				C79701										
				B57801										
				AVSTC-N										
				AVSTC-12										
				PCR-0208-C										
				PCR-2CP-RT-C										